

PROGRAM

**BORT-Quality
Made in Germany**

FOR BORT PEOPLE ARE IN THE CENTER OF INTEREST

For active participation in life – BORT develops products,
oriented acc. to the patients' requirements.

To accelerate recovery and achieve physical harmony.



Only the highest quality for your patients:

The BORT quality management system has been certified acc. to ISO 9001 and EN ISO 13485.

QUALITY MANAGEMENT

certified acc. to
medical products guidelines
EN ISO 13485





QUALITY MANAGEMENT

we are certified
regular voluntary
control acc. to ISO 9001:2000




All products also shown in internet: www.bort.com

- ▶ The easy "product finder" leads you quickly without detours to your requested product!
- ▶ Under "infothek" you will find clinical studies and useful hints on our products.



BORT OFFERS MORE IN MANY INSTANCES

The medical bandages and orthopaedic aids are being developed and tested in close cooperation with physicians and orthopaedic technicians. During the last years this way a very versatile program for the most different applications has been set up. All products are designed for a long-term use and are convincing by an excellent ratio between price and performance.



per palm
Inno app. lomo nel. l'acqua. quanto. e. la. sua. m.
Inno. cimmo. de. capogly. aff. de. fono. de. l'omo. e. il. l'omo. de. l'acqua. e.
m. de. l'acqua. h. l'acqua. de. l'acqua. de. l'acqua. de. l'acqua. de. l'acqua. de. l'acqua.
de. l'acqua. de. l'acqua. de. l'acqua. de. l'acqua. de. l'acqua. de. l'acqua. de. l'acqua.



Neck	6
------	---

Arm	8
-----	---




Finger	16
--------	----

Thumb	20
-------	----



Wrist	26
-------	----



Body	42
------	----

Hernia	52
--------	----

Reha-Products
OTC-Products

Thigh	56
-------	----

Knee	56
------	----

Calf	70
------	----

Ankle	70
-------	----

Foot	75
------	----

Insoles	78
---------	----

Anti-Decubitus	82
----------------	----

Reachers	86
----------	----

Active Supports	88
-----------------	----

ClimaCare [®]	90
------------------------	----

PediSoft [®] /PediSoft [®] TexLine	92
--	----


SoftSocks ergo	94
----------------	----

StabiloColor [®]	95
---------------------------	----

Blister Pad/EasyLife [®]	96
-----------------------------------	----


NewLook	98
---------	----

Foot Rollers	99
--------------	----



No. 127 260
BORT
Cervical Support


Soft open-pore foam, density 45, to wear day and night. Anatomic shape. With recess for larynx. Skin-friendly cover (100 % cotton). Adjustable with Velcro closure. Additional **children's sizes**. Colors: white (No. 127 260), beige (No. 127 360).



No. 127 160
BORT Stabilo® Orthopaedic
Cervical Support

Stable open-pore foam, density 60, to wear by day. With recess for larynx. Adjustable with Velcro closure. Optimal positioning due to anatomic fit. Skin-friendly cover (100 % cotton). Color: beige.

Measuring points:
Eyes looking straight
A = circumference of neck
K = height of chin
(tip of chin to sternoclavicular joint level)



INDICATIONS

For catapult trauma of the cervical vertebra column, so-called whip-lash injury, aching tilted neck, after operations of the cervical vertebra column and after operations of the tilted neck. **Cervical vertebra column syndroms.**

Sizes	Circumference of neck
small	up to 38 cm
medium	38 – 43 cm
large	43 – 48 cm
Height of chin	7,5/9,0/11,0 cm
Children's sizes	Circumf. neck
x-small	up to 33 cm
Height of chin	5,0/6,0 cm


1 cm = 0,3937 inch

INDICATIONS

For operations of tilted neck, after **distortions**, fractures, cervical syndrome as well as catapult trauma and after operations on the cervical vertebra column.

Sizes	Circumference of neck
small	up to 38 cm
medium	38 – 43 cm
large	43 – 48 cm
Height of chin	7,5/9,0/11,0 cm

1 cm = 0,3937 inch



No. 127 560
BORT StabiloNec®
with Reinforcement Brace

Soft, open-pore foam, density 45. With integrated reinforcement splint. Anatomic shape. With recess for larynx. Adjustable with Velcro closure. Cotton cover. Color: white.



No. 127 060
BORT Adjustable
Cervical Collar

Especially good stabilization by non-elastic plastic shell. Adjustable with Velcro closure in height and width. With chin pad (exchangeable). Color: white.


INDICATIONS

Cervical syndrome (grade 2), tilted neck (rheumatic), pains conditions after traumata, distinct catapult trauma, osteochondrosis/spondyloarthrosis of cervical vertebra column, acute radiating neck-shoulder-arm-pain.

Sizes	Circumference of neck
small	up to 38 cm
medium	38 – 43 cm
large	43 – 48 cm

Height of chin 7,5/9,0/11,0 cm
1 cm = 0,3937 inch

Measuring points:
Eyes looking straight
A = circumference of neck
K = height of chin
(tip of chin to sternoclavicular joint level)



INDICATIONS

For injuries and inflammations of the cervical vertebra column, cervical-shoulder-arm-syndrom, painful stress of cervical and shoulder muscles, as well as for rheumatism, arthritis, cervical spondylosis and **catapult trauma**.

Sizes	Circumference of neck
small	up to 40 cm
medium	40 – 44 cm
large	44 – 48 cm

1 cm = 0,3937 inch




No. 121 200
BORT OmoStabil®

Shoulder-arm-orthosis, made of elastic, perforated velours compound material. Cap-shaped recess at shoulder area. Circular with Velcro closure in body area. Free mobility of hand by opening with closing system.

Elastic Velcro closure strap for individual immobilization or positioning of upper arm or shoulder.

Color: white, blue.



Measuring point:
Circumference waist
above crista iliaca.

INDICATIONS

Immobilization/fixing of shoulder and of arm, after contusion of shoulder, after luxation, for periartropathia humero-scapularis (PHS), for injuries of rotary sleeve tendopathy, destruction of acromioclavicular joint, post-operative for immobilization. Subcapitale humerus fractures, after fractures of articular cavity, after fracture of scapule, after arthroscopic operations.

Sizes	Circumference of waist
Size 0	up to 75 cm
Size 1	75 – 95 cm
Size 2	95 – 125 cm
Size 3	125 – 145 cm

Same for left and right.

1 cm = 0,3937 inch



No. 121 100
BORT OmoFX Shoulder-Arm
Abduction Orthosis

Replaceable shoulder-arm abduction orthosis with skin-friendly velours cover, stabilization straps for immobilization of shoulder joint and elbow.

Abduction: approx. 30°.

Optional wedge for increase of abduction to 50° is available (No. 121 160)

Color: blue.

INDICATIONS

Immobilization of shoulder joint and elbow at periarticular irritations and affections in the area of the glenohumeral joint. Calcificating rotary sleeve tendopathy. Bursitis acromialis, supraspinatus syndrome, tendinitis of the long biceps ligament, condition after shoulder corner joint destruction (Tossy I + II).

Sizes	Fits height
Size 1	up to 165 cm
Size 2	over 165 cm

Please indicate left or right.

1 cm = 0,3937 inch



No. 121 900
BORT OmoARS

Arm support cushion for immobilization in 15° or 30° external rotation. Modular system: for resting in 30° external rotation an optional external rotation wedge is available (No. 121 960).

Strapping system ensures optimal fit; simple application by proven closure system. With integrated training ball for autonomous forearm exercises. Color: blue.

INDICATIONS

Post-operative for re-establishment for rotator collar, after subacromial defilee dilation (acc. to Neer), relaxation of the limbus and scapula sections, post-operative immobilization after shoulder prosthesis implantation, conservative immobilization after shoulder luxation, at impingement syndrome and acute and chronic bursitis.

Sizes Fits height

Size 1 up to 165 cm
Size 2 over 165 cm

Same for left and right.

1 cm = 0,3937 inch



No. 121 800
BORT OmoTwinAir

Immobilization of shoulder joint in abduction position of approx. 50°-60°, additional fixation of upper arm by upper arm pad which can be freely positioned. **2-compartment-system to avoid undesired anteversion of upper arm in the shoulder joint while lying.** Especially low net weight. Comfortable hand rest; possibility for exercising lower arm muscular system. Color: blue.

INDICATIONS

Impingement syndrome, resting in abduction position after shoulder joint injuries, condition after rotary sleeve rupture, immobilization after implantation of shoulder prosthesis.

One size only.

Same for left and right.



No. 121 300

BORT OsoTract® Upper Arm-Shoulder-Orthosis

- ▶ Special upper arm brace with 4 weights.
- ▶ 250 g each, to be placed individually into pockets, with Velcro closure, this allows individual weight positioning.
- ▶ Recess in flexion area, therefore comfortable to wear.
- ▶ Anatomic shape provides non-slip position.
- ▶ Color: blue.

INDICATIONS

Impingement syndrome, periarthropathia humero-scapularis (PHS).

One size only.

Same for left and right.




No. 055 100

BORT KubiTal® Elbow Padding Bandage

- ▶ Strong dual-tension quality with special silicone padding.
- ▶ At the elbow the padding is lower, at the edges higher shaped.
- ▶ Padding provides large size protection for the entire elbow area.
- ▶ Color: blue.

*Measuring point:
Circumference of
forearm 5 cm below
elbow joint.*



INDICATIONS

Compression of nervus ulnaris at elbow (cubital tunnel syndrome, sulcus ulnaris), bursitis, after traumatization of olecranon, after surgeries.

Sizes	Measure 5 cm below elbow joint
x-small	22 – 24 cm
small	24 – 26 cm
medium	26 – 28 cm
large	28 – 30 cm
x-large	30 – 32 cm

1 cm = 0,3937 inch



No. 055 300
BORT KubiFX
Elbow-Orthosis

Non-elastic material, inside with terry cloth lining. **With 2 removable alu splints, individually adjustable from 0° to 90°.**

With 4 return straps for solid fit.

Freely positionable elbow pad.

Length: approx. 37 cm.

Color: black.



No. 055 400
BORT KubiFX
long

Non-elastic material, inside with terry cloth lining. **With 2 removable alu splints, individually adjustable from 0° to 90°.** Fixation by means of 5 reversing straps, with Velcro closure.

Elbow padding to be freely positioned.

Length: approx. 53 cm.

Color: black.

INDICATIONS

Immobilization of elbow joint after traumatization and surgeries, bursitis, descending tendopathies, e.g. epicondylitis humeri. Condition after bursectomy.

One size only.

Same for left and right.

INDICATIONS

Immobilization of elbow joint and neutralization of supination and pronation of forearm after traumatization and surgeries, for descending tendopathies, epicondylitis humeri, sulcus ulnaris.

Sizes


Length Olecranon to metacarpal knuckle point V


Size 1	28 – 31 cm
Size 2	31 – 34 cm
Size 3	34 – 38 cm

Please indicate left or right.

1 cm = 0,3937 inch

*Measuring point:
Olecranon to
metacarpal
knuckle point V.*






No. 054 200
BORT Elastic
Elbow Support

Circular knitted, thus without seams. Dual-tension elastic, does not slip, thin material, slip-on. Compression class II. Color: beige.

 **Microfibre comfort**




No. 104 200
BORT
Elbow Bandage

Strong but nevertheless comfortable bandage with good wearing comfort. Anatomic shape-knitted for optimised fit. Color: beige.

No. 055 200
BORT Dual-Tension
Elbow Bandage

With one stay on each side. Color: silvergrey.



Measuring point:
A = Circumference of
elbow joint.

INDICATIONS

Articular effusions, after operations at elbow as well as for **sports injuries**. Distortions, contusions, arthrosis.

Sizes	Circumference of elbow
x-small (No. 054 200)	22 – 24 cm
small	24 – 26 cm
medium	26 – 28 cm
large	28 – 30 cm
x-large	30 – 32 cm

1 cm = 0,3937 inch

INDICATIONS

Articular effusion, after operations at elbow as well as **sports injuries**, distortions, contusions, arthrosis.

Sizes	Circumference of elbow
small	up to 26 cm
medium	26 – 28 cm
large	28 – 30 cm

1 cm = 0,3937 inch



No. 122 500

BORT EpiContur® Epicondylitis-Brace with 2 pads


- ▶ Epicondylitis-brace with specific fit for friction or for compression of ligament insertions.
- ▶ Conic shape, therefore fixation close to joint possible.
- ▶ With 2 freely positionable silicone pads.
- ▶ Individually adjustable pressure by reversed strap with Velcro closure.
- ▶ Colors: blue, beige.

No. 022 500

BORT EpiContur® with 1 pad

Model as above, however only with 1 silicone pad.
Colors: blue, beige.

*Measuring point:
Circumference of
forearm 5 cm below
elbow joint.*

**INDICATIONS**

Epicondylitis. Against aching elbows as a sponge in the case of sports injuries, overstrain, tennis elbow, golf elbow.

Sizes**Measure
5 cm below
elbow joint**

Size 0	21 – 25 cm
Size 1	25 – 30 cm
Size 2	30 – 35 cm

1 cm = 0,3937 inch



No. 122 580
**BORT Epicondylitis-Brace
 with Ulnar Relief**

With 2 pads, aluminium strap for defined pressure action without undesired compression, specially formed to avoid pressure to the ulna. Color: white.



No. 122 600
**BORT Epicondylitis Bandage
 with Padding**

Strong dual-tension quality. With two structured silicone pads for intermittent massage. Shape-knitted with pressure diminishing zone above olecranon. Free mobility of elbow joint. Colors: beige, silvergrey.

 **Microfibre comfort**

INDICATIONS

Epicondylitis. For defined pressure action without participation of ulna.

Sizes	Measure 5 cm below elbow joint
Size 0	23 – 25 cm
Size 1	25 – 27 cm
Size 2	27 – 29 cm


1 cm = 0,3937 inch

INDICATIONS


Epicondylitis (acute or chronic), joint arthrosis, condition after fractures near the joint, articular effusion.

Sizes	Measure 5 cm below elbow joint
x-small	22 – 24 cm
small	24 – 26 cm
medium	26 – 28 cm
large	28 – 30 cm
x-large	30 – 32 cm

1 cm = 0,3937 inch



*Measuring point:
 Circumference of
 forearm 5 cm below
 elbow joint.*



No. 112 560

BORT MobiDig® Flexion Finger Splint

- ▶ Finger splint for flexion of DIP or PIP joints.
- ▶ High wearing comfort by ergonomic, wear-resistant padding.
- ▶ By careful bending of spring, flexion can be increased or decreased.
- ▶ Linear spring force.


Measuring points:

Flexion **PIP joint**:

L_1 = Finger basis joint to distal finger fold.

Flexion **DIP joint**:

L_2 = PIP joint to distale finger tip.

**INDICATIONS**


Flexion of PIP(L_1)- or DIP(L_2)-joints.

Sizes**Length**

L_1 = Flexion PIP-joint
or L_2 = Flexion DIP-joint

x-small	up to 3,3 cm
small	3,3 – 4,6 cm
medium	4,6 – 5,2 cm
large	5,2 – 6,0 cm
x-large	over 6,0 cm

1 cm = 0,3937 inch



No. 112 510

BORT MobiDig® Extension Finger Splint

- ▶ Finger splint for extension of PIP joints.
- ▶ High wearing comfort by ergonomic, wear-resistant padding.
- ▶ By carefully bending spring, tension can be increased or decreased.
- ▶ Proximal padding positioned in area of palm of hand.
- ▶ Linear spring force.
- ▶ Color: black.



Measuring point:
L = Distal fold of palm of hand
to distal finger fold.



Measuring point:
Finger width at
PIP joint

INDICATIONS

Extension of PIP joint.


Sizes	Length (Distal fold of palm to distal finger fold)	Finger width at PIP joint
small	6,3 – 7,6 cm	1,5 – 1,7 cm
small plus	6,3 – 7,6 cm	1,7 – 1,9 cm
medium	7,6 – 8,9 cm	1,7 – 1,9 cm
medium plus	7,6 – 8,9 cm	1,9 – 2,1 cm
large	8,9 – 10,2 cm	1,9 – 2,1 cm
large plus	8,9 – 10,2 cm	2,1 – 2,3 cm

1 cm = 0,3937 inch



No. 112 520
BORT Capener Splint
Extension

For extension of PIP joints. With lateral spiral springs for anatomically correct movement. Soft padded. With Velcro closure for fixing. Progressive spring force.



No. 112 590
BORT Capener Splint
Flexion

With two lateral spiral springs for anatomically correct movement. Large dorsal pad to spread pressure. Proximal strap for appropriate tension. Progressive spring force.

Measuring point:
L = Distal fold of palm
of hand to distal
finger fold.



INDICATIONS

Extension of PIP joints.

Sizes

Length

(Distal fold of palm of hand
to distal finger fold)

small	6,3 – 7,0 cm
medium	7,0 – 7,6 cm
large	7,6 – 8,0 cm

1 cm = 0,3937 inch

Measuring point:
L = Metacarpal
knuckle to distal
finger fold.



INDICATIONS

Flexion of PIP joints.

Sizes

Length

(Metacarpal knuckle
to distal finger fold)

small	6,3 – 6,9 cm
medium	6,9 – 7,5 cm
large	7,5 – 8,2 cm
x-large	8,2 – 8,5 cm

1 cm = 0,3937 inch



No. 112 770
BORT DigiSoft® - Finger Orthosis

- ▶ Finger splint for fixation of 1 or 2 fingers.
- ▶ Pad with soft, ventilating material.
- ▶ Alu-splint can be remodelled.
- ▶ 3 Velcro closures for optimal fixation, center Velcro closure can be placed in variable position.
- ▶ Color: blue.

Measuring point:
L = Begin of
palm to end
of finger.



INDICATIONS

After injuries and operations, immobilization after plaster cast. Redressement.


Sizes	Length of hand
Size 1	up to 14 cm
Size 2	over 14 cm

1 cm = 0,3937 inch



No. 112 710
BORT SOFT
Thumb Splint, short

Comfortable soft thumb splint. Dry wearing feeling due to naturally breathing and humidity draining material. Integrated anatomically pre-shaped splint for partial immobilization of thumb. Free mobility of the fingers. Color: blue/black.



No. 112 720
BORT SOFT
Thumb Splint, long

Comfortable soft thumb splint. Dry wearing feeling due to naturally breathing and humidity draining material. Integrated anatomically pre-shaped splint for partial immobilization of thumb. Free mobility of the fingers. Extended splint up to wrist area. Color: blue/black.

INDICATIONS

After thumb fractures, arthritis, arthrosis of thumb saddle joint and/or metacarpal phalangeal joint, ulnar ligament lesions (ski thumb).

Sizes	Circumference of wrist
small	15 – 17 cm
medium.	17 – 19 cm
large.	19 – 21 cm

Same for left and right.

1 cm = 0,3937 inch

INDICATIONS

After thumb fractures, arthritis, arthrosis of thumb saddle joint and/or metacarpal phalangeal joint, ulnar/radial ligament injury of thumb. Prevents false movements, late stage after sinew reconstruction, tendopathy.

Sizes	Circumference of wrist
small	15 – 17 cm
medium.	17 – 19 cm
large.	19 – 21 cm

Same for left and right.

1 cm = 0,3937 inch

Measuring point:
H = Circumference
of wrist proximal
of styloideus.





No. 112 420
BORT
 Thumb Wrist Bandage

Anatomically form-knitted. Compression class II. Supports wrist and thumb. 2 lateral stays and 1 adjustable dorsal splint in thumb area. With additional Stabilo® strap for better support.
 Colors: beige, blue.



No. 112 050
BORT
 SellaTex®

Rigid bandage with embraced thumb. Comfortable and dry wearing impression by naturally breathing and humidity discharging material. With three Velcro closures and volar wrist splint. Dorsal aluminium thumb splint with additional strong, circular Velcro strap, approx. 20 cm long. Shorter volar splint in children's size. Color: beige, black/blue.

INDICATIONS

After thumb fractures, arthritis, ulnar/radial ligament injury of thumb.

Measuring point:
H = Circumference of wrist proximal of styloideus.



INDICATIONS

After thumb fracture, arthritis, ulnar/radial ligament injury of thumb. Protection against wrong movements, late state after ligament reconstruction, tendopathy, arthrosis of thumb saddle joint and/or metacarpal phalangeal joint.

Sizes	Circumference of wrist
small	up to 15 cm
medium	15 – 17 cm
large	17 – 19 cm
x-large	19 – 21 cm

Same for left and right.

1 cm = 0,3937 inch

Sizes	Circumference of wrist
For children (blue only)	
x-small	13 – 15 cm
small	15 – 17 cm
medium	17 – 19 cm
large	19 – 21 cm
x-large	21 – 23 cm

Please indicate left or right.

1 cm = 0,3937 inch



No. 112 080
BORT SellaFix® K
Thumbstall

Short version in the carpal joint area for free mobility of wrist. Can be reworked individually with hair-dryer and scissors. Perforation provides for venting. Adjustable by Velcro closure. Terry padding strips can be positioned individually. Inner face completely flocked with textile. Color: beige.



No. 112 900
BORT
SellaFix® P

Plus extended thumb tube and Velcro closure for immobilization of thumb joint, can be readjusted individually with hair-dryer and scissors. Perforation provides for venting. Adjustable by circular Velcro closure, with terry padding strips for individual positioning. Inner face completely flocked with textile. Color: beige.

INDICATIONS

For immobilization of - metacarpal phalangeal joint
with full mobility of wrist, in case of strains and sprains.

Sizes	Circumference of wrist
x-small	up to 15 cm
small	15 – 17 cm
medium	17 – 19 cm
large	19 – 21 cm
x-large	21 – 23 cm

Please indicate left or right.

1 cm = 0,3937 inch

INDICATIONS

For immobilization of - metacarpal phalangeal joint - thumb saddle joint - interphalangeal joint
Tendopathies, fractures.

Sizes	Circumference of wrist
x-small	up to 15 cm
small	15 – 17 cm
medium	17 – 19 cm
large	19 – 21 cm
x-large	21 – 23 cm

Please indicate left or right.

1 cm = 0,3937 inch

Measuring point:
H = Circumference
of wrist
proximal of
styloideus.





No. 112 040
BORT SellaFix[®] N Metacarpus Splint

- ▶ Thumb saddle joint orthosis made of thermoplastic material.
- ▶ Easy correction with hair-dryer and scissors.
- ▶ This allows individual fitting.
- ▶ Fine venting perforation.
- ▶ Washable.
- ▶ Terry padding can be individually positioned.
- ▶ Inner face completely flocked with textile.
- ▶ Color: beige.

Measuring point:
H = Circumference of wrist proximal of styloideus.



INDICATIONS

For immobilization of
 - metacarpal phalangeal joint and
 - thumb saddle joint

Ulnar ligament lesion (ski thumb), partial immobilization, rhizo-arthritis.

Sizes	Circumference of wrist
x-small	up to 15 cm
small	15 – 17 cm
medium	17 – 19 cm
large	19 – 21 cm
x-large	21 – 23 cm

Please indicate left or right.

1 cm = 0,3937 inch



No. 105 600
BORT SellaDur® FLOC

- ▶ Thumb saddle joint orthosis with excellent fit.
- ▶ In the critical thumb-forefinger-area special high formed edge to avoid pressure spots.
- ▶ Additional **flocking** with textiles provides optimal wearing comfort.
- ▶ Perforation provides for venting.
- ▶ Washable.
- ▶ Color: grey.

Measuring point:
H = Circumference
of wrist
proximal of
styloideus.



INDICATIONS

For immobilization of
- metacarpal phalangeal joint and
- thumb saddle joint

Ulnar ligament lesions (ski thumb), partial immobilization, rhizo-arthritis.

Sizes	Circumference of wrist
small	15 – 17 cm
medium.....	17 – 19 cm
large.....	19 – 21 cm

Please indicate left or right.

1 cm = 0,3937 inch



No. 105 400
BORT
SellaDur®

Thumb saddle joint orthosis made of special tough Duroplast. With perforation, therefore very comfortable to wear. In the critical thumb-forefinger-area special high formed edge to avoid pressure spots. Pads to be individually positioned. Color: blue.



No. 103 500
BORT Thumb/
Hand Orthosis

Made of thermoplastic material. Easy correction with hair-dryer at 60°C, allowing individual adjustment. Fine venting perforation. Inner face completely flocked with textile. Washable. Color: beige.

INDICATIONS

For immobilization of - metacarpal phalangeal joint and - thumb saddle joint

Ulnar ligament lesion (ski thumb), partial immobilization, rhizoarthrosis.

Sizes	Circumference of wrist
small	15 – 17 cm
medium.	17 – 19 cm
large.	19 – 21 cm

Please indicate left or right.

1 cm = 0,3937 inch

INDICATIONS

For immobilization of metacarpal phalangeal joint, thumb saddle joint and interphalangeal joint.

Immobilization of wrist. Strains and sprains.

Sizes	Circumference of wrist
x-small	up to 15 cm
small	15 – 17 cm
medium.	17 – 19 cm
large.	19 – 21 cm
x-large	21 – 23 cm

Please indicate left or right.

1 cm = 0,3937 inch

Measuring point:
H = Circumference
of wrist
proximal of
styloideus.





No. 112 130
BORT
Wrist Strap

Made of leather, with thumb loop and circumferential forked strap, width 6 cm, adjustable, inner face lined.
Color: beige.



No. 112 110
BORT
Wrist Bandage

Elastic bandage to support wrist for easy pains. With Velcro closure, width 8 cm.
Colors: beige, black, multicolored.

INDICATIONS

Distortions and **overstressing of wrist joint.**

Sizes	Circumference of wrist
xx-small	15 cm
x-small	16 cm
small	17 cm
medium	18 cm
large	19 cm
x-large	20 cm

Please indicate left or right.

1 cm = 0,3937 inch


INDICATIONS

For tendosynovitis at wrist, strains of wrist and after removal of plaster.

Sizes	Circumference of wrist
Size 1	up to 17 cm
Size 2	17 – 19 cm
Size 3	19 – 22 cm

1 cm = 0,3937 inch

Measuring point:
H = Circumference of wrist proximal of styloideus.





No. 104 300
BORT Wrist Support
Closed Form

Circular wrist support without closing, slip-on version. Ideal for early overstress symptoms. Color: beige.



No. 112 010
BORT Stabilo®
Wrist Bandage

Elastic bandage with Velcro closure and additional Stabilo® strap for increased support effect. Individually adjustable, width 8 cm. Colors: beige, white, black, white with blue strap.

Wrist

INDICATIONS

For tendosynovitis at wrist, strains of wrist and after removal of plaster.

Sizes	Circumference of wrist
small	up to 17 cm
medium.	17 – 19 cm
large.	19 – 22 cm

1 cm = 0,3937 inch


INDICATIONS

For tendosynovitis, strains of the wrist and **after removal of plaster.**

Sizes	Circumference of wrist
Size 1	up to 17 cm
Size 2	17 – 19 cm
Size 3	19 – 22 cm

1 cm = 0,3937 inch

Measuring point:
H = Circumference of wrist proximal of styloideus.





No. 182 520
BORT StabiloPren
Wrist Bandage

Light neoprene bandage with Velcro closure, thumb loop avoids displacing.
Color: black.



No. 112 020
BORT Wrist Support
with edged Thumb opening

Optimized fit in palmar area. Anatomically form-knitted with edged thumb opening, therefore non-slipping fit.
Compression class II.
Color: beige.

Measuring point:
 H = Circumference
of wrist
proximal of
styloideus.



INDICATIONS

After distortions and **overstressing of wrist joint.**

One size only.

Same for left and right.

INDICATIONS

For luxations and minor distortions of the wrist. **After removal of plaster** or immobilizing treatment, tendovaginitis.

Sizes

Circumference of wrist

x-small	up to 13 cm
small	13 – 15 cm
medium	15 – 17 cm
large	17 – 19 cm
x-large	19 – 21 cm

Same for left and right.

1 cm = 0,3937 inch



No. 112 120

Wrist Support with Thumb opening, additional Strap

- ▶ Optimized fit in palmar area.
- ▶ Anatomically form-knitted.
- ▶ Supports wrist up to back of the hand.
- ▶ With edged thumb opening, therefore non-slipping fit.
- ▶ Additional Stabilo® strap for effective support.
- ▶ Color: beige.

Measuring point:
H = Circumference
 of wrist
 proximal of
 styloideus.



INDIKATIONEN

For luxations and minor distortions. Strains of the wrist, **after removal of plaster**, tendovaginitis.

Sizes	Circumference of wrist
x-small	up to 13 cm
small	13 – 15 cm
medium	15 – 17 cm
large	17 – 19 cm
x-large	19 – 21 cm

Please indicate left or right.

1 cm = 0,3937 inch



No. 112 220

BORT ManuZip® volar Medianus Splint

- ▶ Anatomically form-knitted.
- ▶ **Volar** aluminium splint for individual adjustment.
- ▶ Viscoelastic silicone pad with intermitting massage effect.
- ▶ **Extended zipper up to thumb, therefore very easy to apply.**
- ▶ With integrated slip-on aid for easy putting on, also if gripping ability is restricted.
- ▶ Especially large zipper carriage for easy gripping with locking device; ideal for rheumatics.
- ▶ Additional Stabilo® strap for safe stabilization.
- ▶ Colors: beige, blue.

Measuring point:
*H = Circumference
 of wrist
 proximal of
 styloideus.*



INDICATIONS

For arthrosis of carpus, distortions, luxations and functionally conditioned immobilization. Irritations, arthrosis, scaphoid-pseudarthrosis, after removal of plaster, arthralgias, rheumatism and drop hand.

Sizes

Circumference of wrist

x-small	up to 13 cm
small	13 – 15 cm
medium	15 – 17 cm
large	17 – 19 cm
x-large	19 – 21 cm

Please indicate left or right.

1 cm = 0,3937 inch



No. 112 620
BORT ManuZip® Eco

Anatomically form-knitted. **Volar** aluminium splint for individual adjustment. **Extended zip closure up to thumb**, therefore easy to put on. Additional Stabilo® strap for secure stabilization. Colors: beige, blue



No. 112 320
BORT ManuZip® dorsal

Anatomically form-knitted. **Dorsal** aluminium splint for individual adjustment. 2 viscoelastic silicone pads with intermitting massage effect. **Extended zipper up to thumb, therefore very easy to apply.**

With integrated slip-on aid for easy putting on, also if gripping ability is restricted. Especially large zipper carriage for easy gripping with locking device; ideal for rheumatics. Additional Stabilo® straps for safe stabilization. Color: beige.

Measuring point:
H = Circumference of wrist proximal of styloideus.



INDICATIONS

For arthrosis of carpus, distortions, luxations and functionally conditioned immobilization. Irritations, arthrosis, scaphoid-pseudarthroses. After removal of plaster, arthralgies, rheumatism, drop hand.

INDICATIONS

For arthrosis of carpus, distortions, luxations and functionally conditioned immobilization, irritations, arthrosis, scaphoid-pseudarthrosis. After removal of plaster, arthralgies, rheumatism.

Sizes	Circumference of wrist
x-small	up to 13 cm
small	13 – 15 cm
medium.	15 – 17 cm
large.	17 – 19 cm
x-large	19 – 21 cm

Please indicate left or right.

1 cm = 0,3937 inch

Sizes	Circumference of wrist
x-small	up to 13 cm
small	13 – 15 cm
medium.	15 – 17 cm
large.	17 – 19 cm
x-large	19 – 21 cm

Please indicate left or right.

1 cm = 0,3937 inch



No. 103 300

BORT Wrist Brace with Alu-Splint

- ▶ Completely to open, with Velcro closures at back of the hand.
- ▶ With integrated slip-on aid for easy putting on, also if gripping ability is restricted.
- ▶ With moldable aluminium splint, 16 cm long.
- ▶ Splint in durable nylon jacket.
- ▶ Colors: beige, blue/black.

No. 103 301

BORT Wrist Brace with additional strap

Model as above, however with additional strap around wrist.
Colors: beige, blue/black.

Measuring point:
 H = Circumference
of wrist
proximal of
styloideus.

**INDICATIONS**

Carpal tunnel syndrome, tendosynovitis, rheumatism, strains of the wrist and after removal of plaster, immobilization of the wrist after non-bone and affected injuries, radial nerve paralysis (drop hand).

Sizes**Circumference of wrist**

x-small	up to 15 cm
small	15 – 17 cm
medium	17 – 19 cm
large	19 – 21 cm
x-large	21 – 23 cm

Please indicate left or right.

1 cm = 0,3937 inch



No. 103 360
BORT Wrist
and Arm Brace

Volar forearm splint with long, moldable aluminium splint, **22 cm long**, brace can be opened completely. With Velcro closures at back of the hand. With integrated slip-on aid for easy putting on, also if gripping ability is restricted.
Colors: beige, blue/black.



No. 103 580
BORT Arm- and
Wrist Splint

Volar forearm splint with strap and long, moldable aluminium splint, **22 cm long**, brace can be opened completely. With Velcro closures at back of the hand. **With additional Stabilo® Strap.** With integrated slip-on aid for easy putting on, also if gripping ability is restricted.
Colors: beige, blue/black.

INDICATIONS

Carpal tunnel syndrome, tendosynovitis, rheumatism, strains of the wrist and after removal of plaster, immobilization of the wrist after non-bone and affected injuries, radial nerve paralysis.

Sizes	Circumference of wrist
x-small	up to 15 cm
small	15 – 17 cm
medium	17 – 19 cm
large	19 – 21 cm
x-large	21 – 23 cm


Please indicate left or right.

1 cm = 0,3937 inch

INDICATIONS

Carpal tunnel syndrome, tendosynovitis, strains of the wrist and after removal of plaster, immobilization of the wrist after non-bone and affected injuries, radial nerve paralysis, stabilization of carpus.

Measuring point:
H = Circumference of wrist proximal of styloideus.



Sizes	Circumference of wrist
x-small	up to 15 cm
small	15 – 17 cm
medium	17 – 19 cm
large	19 – 21 cm
x-large	21 – 23 cm

Please indicate left or right.

1 cm = 0,3937 inch



No. 112 730
BORT
ManuStabil® short

Solid wrist bandage. Comfortable and dry wearing impression by naturally breathing and humidity discharging material. With three Velcro closures and volar and dorsal splint and additional strong circular Velcro strap.

Approx. **20 cm** long.

Shorter volar and dorsal splint in children's size.

Colors: beige, black/blue.



No. 112 750
BORT SOFT Volar
CTS Splint

Padded with soft naturally breathing material. Comfortable to wear. Moldable volar splint. 3 soft padded Velcro closure straps for safe fit. Free movability of fingers.

Splint length 20 cm.

Color: blue.

Measuring point:
 H = Circumference
of wrist
proximal of
styloideus.



INDICATIONS

Carpal tunnel syndrome, tendosynovitis, rheumatism, strains of wrist and after removal of plaster. Immobilization of wrist after non-bone affected injuries, radial nerve paralysis (drop hand). **Also to wear at night.**

Sizes	Circumference of wrist
-------	------------------------

For children (blue only)	13 – 15 cm
x-small	up to 15 cm
small	15 – 17 cm
medium	17 – 19 cm
large	19 – 21 cm
x-large	21 – 23 cm

Please indicate left or right.

1 cm = 0,3937 inch

INDICATIONS

Carpal tunnel syndrome, tendon injuries, wrist sprains, after fractures.

Sizes	Circumference of wrist
-------	------------------------

small	15 – 17 cm
medium	17 – 19 cm
large	19 – 21 cm


Please indicate left or right.

1 cm = 0,3937 inch



No. 182 500
BORT StabiloPren
 Arm and Wrist Support

Made of Neopren. Optimized fit by anatomic shape. Inside with jersey cloth laminated, therefore skin-friendly. With integrated slip-on assist and reversing rein for easy application, even with restricted gripping ability. **25 cm** long. Splint moldable.
 Color: black/blue.



No. 102 200
BORT ManuStabil®

Solid wrist orthosis with volar and dorsal splint, **25 cm** long. Easy and quick closing by special straps, secured by Velcro closure.
 Color: blue.

INDICATIONS

Strains, sprains, carpal tunnel syndrome.

Sizes	Circumference of wrist
x-small	up to 16 cm
small	16 – 18 cm
medium.	18 – 20 cm
large.	20 – 22 cm
x-large	22 – 24 cm

Please indicate left or right.

1 cm = 0,3937 inch

INDICATIONS


Carpal tunnel syndrome, tendosynovitis, strains of the wrist and after removal of plaster. For immobilization of the wrist after bone and non-bone injuries, radial nerve paralysis, stabilization of carpus.

Sizes	Circumference of wrist
x-small	up to 15 cm
small	15 – 17 cm
medium.	17 – 19 cm
large.	19 – 24 cm

Please indicate left or right.

1 cm = 0,3937 inch

Measuring point:
H = Circumference of wrist proximal of styloideus.






No. 102 300

BORT ManuDisc based on Dr. Schütz

- ▶ Firm wrist bandage with volar T-splint, **ulnar embracing and bilateral support**.
- ▶ **Multidirectional stability** in the wrist at free mobility of fingers.
- ▶ **Comfortable and dry-wearing sensation by air ventilating and humidity deducting material.**
- ▶ Soft padding in problem areas (e.g. commissur D1/D2.)
- ▶ Individually adjustable.
- ▶ Color: blue/black.

Measuring point:
H = Circumference
of wrist
proximal of
styloideus.



INDICATIONS

Carpal tunnel syndrome, **disc injuries**, arthrosis of carpus, ligament lesions in carpal area, **post-operative immobilization**, distortions, arthralgies, rheumatism, drop hand, condition after metacarpus fracture for functional therapy.

Sizes

Circumference of wrist

small	15 – 17 cm
medium	17 – 19 cm
large	19 – 21 cm
x-large	21 – 23 cm

Please indicate left or right.

1 cm = 0,3937 inch



No. 103 200
BORT ManuCarpal®-Orthosis

Very light (approx. 100 g), unobstrusive splint. Full mobility of the fingers. No undesired compression, with padding. Inner face completely flocked with textile, therefore comfortable to wear. Made of low temperature thermoplastic material. Adjustment easily possible with hair-dryer and scissors. Therefore individual readjustment and corrections possible in case of anatomic changes. Washable. Color: beige.



No. 102 100
BORT ManuCarpal®-SOFT
Orthosis

Individually mouldable orthosis, padded with soft naturally breathing material. Wide finger support, simple fixture of thumb. 4 Velcro closures for good fixation. Color: blue.

Measuring point:
H = Circumference
of wrist
proximal of
styloideus.



INDIKATIONEN

Carpal tunnel syndrome. Tendopathies, immobilization of wrist.

Sizes	Circumference of wrist
x-small	up to 15 cm
small	15 – 17 cm
medium	17 – 19 cm
large	19 – 21 cm
x-large	21 – 23 cm

Please indicate left or right.

1 cm = 0,3937 inch

INDIKATIONEN

For night resting, carpal tunnel syndrome, immobilization of wrist after scaphoid fractures, resting in case of tendosynovitis, lymphangitis and post-operative.

Sizes	Circumference of wrist
small	15 – 17 cm
medium	17 – 19 cm
large	19 – 21 cm

Please indicate left or right.

1 cm = 0,3937 inch



No. 133 300
BORT
Finger Joint Support

With long aluminium splint, mouldable, adjustable with Velcro closure. Comfortable and dry wearing impression by naturally breathing and humidity discharging material. With additional Stabilo® strap for safe stabilization.

With integrated slip-on aid for easy putting on, also if gripping ability is restricted.

Colors: beige, blue/black.



No. 112 760
BORT SOFT
Dorsal CTS Splint

With soft, naturally breathing material padded dorsal splint. Fixed with 3 Velcros. Palmar bow supports wrist. **Gripping ability maintained.** Easy to readjust. Color: blue.

INDICATIONS

Carpal tunnel syndrome, tendosynovitis, strains of the wrist and after removal of plaster. For immobilization of the wrist after non-bone affected injuries, radial nerve paralysis, stabilization of carpus. Temporary immobilization in case of hand and finger injuries, tendovaginitis of fingers and hand.

Sizes	Circumference of wrist
x-small	up to 15 cm
small	15 – 17 cm
medium	17 – 19 cm
large	19 – 21 cm
x-large	21 – 23 cm

Please indicate left or right.

1 cm = 0,3937 inch

Measuring point:
H = Circumference
of wrist
proximal of
styloideus.



INDICATIONS

Carpal tunnel syndrome, ligament injuries, wrist sprains, after fractures.

Sizes	Circumference of wrist
small	15 – 17 cm
medium	17 – 19 cm
large	19 – 21 cm

Please indicate left or right.


1 cm = 0,3937 inch



No. 103 800
BORT Dorsal CTS/Radialis Splint

- ▶ Functional positioning of the hand.
- ▶ Full gripping ability.
- ▶ Made of thermoplastic material.
- ▶ Easy adjustment with hair-dryer at 60°C.
- ▶ Therefore individual readjustment and corrections possible in case of anatomic changes.
- ▶ Easily to be put on by patient.
- ▶ Inner face completely flocked with textile.
- ▶ Washable.
- ▶ Color: beige.

Measuring point:
H = Circumference
of wrist
proximal of
styloideus.



INDICATIONS

Radial nerve paralysis (drop hand), **carpal tunnel syndrome.**

Sizes	Circumference of wrist
x-small	up to 15 cm
small	15 – 17 cm
medium	17 – 19 cm
large	19 – 21 cm
x-large	21 – 23 cm

Please indicate left or right.

1 cm = 0,3937 inch



No. 103 700
BORT Paresis Splint

- ▶ Positioning in the intrinsic-plus position.
- ▶ Made of thermoplastic material.
- ▶ Adjustment easily possible with hair-dryer at 60°C.
- ▶ Therefore individual readjustment and corrections possible in case of anatomic changes.
- ▶ Perforated and padded.
- ▶ Inner face completely flocked with textile.
- ▶ Washable.
- ▶ Color: beige.

Measuring point:
H = Circumference
of wrist
proximal of
styloideus.



INDICATIONS

Peripheric nerve paralysis, apoplexy, post-traumatic sinew injury, carpal tunnel syndrome.

Sizes

Circumference of wrist

x-small	up to 15 cm
small	15 – 17 cm
medium	17 – 19 cm
large	19 – 21 cm
x-large	21 – 23 cm

Please indicate left or right.

1 cm = 0,3937 inch



No. 103 100
BORT
Resting Splint

Extra light resting splint for wrist, finger and thumb. Deformable metal splints covered with soft, lasting, foamed plastics. Washable. Can immediately be adjusted individually without tools.
Color: anthracite.



No. 103 900
BORT
Rheumatism Splint

Splint for positioning at night.

Support of the hand vaults with immobilization of the thumb. Made of thermoplastic material. Easy adjustment with hair-dryer. Inner face completely flocked with textile.
Color: beige.

INDICATIONS

Positioning splint for wrist, fingers and thumb.

Sizes	Circumference of wrist
small	15 – 17 cm
medium	17 – 19 cm
large	19 – 21 cm

Please indicate left or right.

1 cm = 0,3937 inch

INDICATIONS

Arthritic deformities – prevention and correction.


Anthalgic immobilization. Support of the hand vault.

Sizes	Circumference of wrist
x-small	up to 15 cm
small	15 – 17 cm
medium	17 – 19 cm
large	19 – 21 cm
x-large	21 – 23 cm

Please indicate left or right.

1 cm = 0,3937 inch

Measuring point:
H = Circumference of wrist proximal of styloideus.





No. 102 900
BORT Rib Belt
for Women


Elastic back part combined with rigid lateral and front parts. Wide Velcro system allows individual adjustment. Comfortable wearing features by cotton satin and fine terry cloth. Height: dorsal 16 cm, ventral 12 cm. Color: white.




No. 102 800
BORT Rib Belt
for Men

Elastic back part combined with rigid lateral and front parts. Wide Velcro system allows individual adjustment. Comfortable wearing features by cotton satin and fine terry cloth. Height: 16 cm. Color: white.

Measuring point:
Circumference
of chest.
B = women.



Measuring point:
Circumference
of chest.
A = men.



INDICATIONS

For contusions and strains of intermediate rib muscles, **rib contusions, rib fractures**, sternum contusions and sternum fractures.

Sizes	Circumference of chest (B)
small	up to 66 cm
medium	.66 – 76 cm
large	.76 – 89 cm
x-large	.89 – 102 cm

1 cm = 0,3937 inch

INDICATIONS

For contusions and strains of intermediate rib muscles, **rib contusions, rib fractures**, sternum contusions and sternum fractures.

Sizes	Circumference of chest (A)
small	up to 82 cm
medium	.82 – 96 cm
large	.96 – 112 cm
x-large	.112 – 128 cm

1 cm = 0,3937 inch




No. 104 600
BORT StabiloFix®

- ▶ Application as support for **up-right position** for posture correction.
- ▶ Application as **clavicle bandage** in case of clavicular injuries.
- ▶ Soft, non-elastic, skin-friendly material for good wearing comfort.
- ▶ With tucks for easy shortening of straps.
- ▶ Adapts optimally to body shape, especially in critical armpit area.
- ▶ Colors: white, blue.

Body
Hernia

Measuring point:
Circumference
of chest.
B = women.



INDICATIONS

For postural disturbance, collarbone fracture.

Sizes	Circumference of chest (B)
For children (blue only)	50 – 70 cm
For adults	over 70 cm

1 cm = 0,3937 inch



No. 112 610


BORT Vario Back Support with Pelotte

- ▶ In back area non-elastic cotton drill.
- ▶ Elastic pull straps from dorsal to ventral.
- ▶ 3 supporting splints in dorsal area.
- ▶ With lumbo sacral friction pelotte.
- ▶ Height: ventral 17 cm, dorsal 24 cm.
- ▶ Colors: beige, black.

No. 112 650

BORT Vario Back Support without Pelotte

Model as above, however without pelotte.
Colors: beige, black.



Measuring point:
Circumference waist
above *Crista iliaca*.


INDICATIONS

Back muscle tensions, spondylolyse, pains in the vertebrae, lumbago, intervertebral disc problems, spondylose, chronic pains in the spine, pains in sacroiliac joint after intervertebral disc surgeries, rheumatic problems of lumbar vertebrae, irritations of nerve roots.

Sizes**Circumference of waist**

small	up to 89 cm
medium	89 – 104 cm
large	104 – 119 cm
x-large	119 – 134 cm

1 cm = 0,3937 inch




No. 104 710

BORT Stabilo® Back Bandage with Pelotte

- ▶ Anatomically form-knitted bandage.
- ▶ With lumbo sacral pelotte made of silicone.
- ▶ 2 anatomical preformed splints in dorsal area.
- ▶ Easy to put on by hand straps.
- ▶ Height: ventral 17 cm, dorsal 26 cm.
- ▶ Colors: beige, silvergrey.

No. 104 700
BORT Stabilo® Back Bandage without Pelotte

Model as above, however with 4 splints and without pelotte.
Colors: beige, silvergrey.



*Measuring point:
Circumference waist
above Crista iliaca.*

INDICATIONS

For back muscle tensions, spondylolyse, spondylolisthesis, pains in the vertebrae, lumbago, damaged discs, spondylose, **chronical pains in the spine**, after intervertebral disc surgeries.


Sizes	Circumference of waist
small	up to 89 cm
medium	.89 – 104 cm
large	.104 – 119 cm
x-large	.119 – 134 cm

1 cm = 0,3937 inch



No. 104 690
BORT Stabilo Basic
Back Support with Pelotte

Active bandage for muscular stabilization of the lumbar vertebral spine. Anatomical formknitted bandage with high wearing comfort. Comfortable and low effort applying of bandage by hand grip loops. 4 splints. With sacral pelotte which can be variably positioned by Velcro fitting. Height: ventral 17 cm, dorsal 22 cm. Colors: beige, silvergrey.




Measuring point:
Circumference waist
above *Crista iliaca*.

INDICATIONS

Dorsalgia, lumbalgia, lumbago, ischi-
algia, irritation of nerve roots.

Sizes	Circumference of waist
small	up to 89 cm
medium	89 – 104 cm
large	104 – 119 cm
x-large	119 – 134 cm

1 cm = 0,3937 inch



No. 104 730
BORT StabiloPlus
Bridging Bandage

Relief of affected vertebral spine section by relief of muscular system and bridging of lumbar lordosis as well as activation of abdominal muscular system with concurrent straightening of the pelvis. The degree of stabilization is adjustable to the course of therapy. The skin-friendly knit-work transfers quickly humidity of the skin surface to the outside and cares for comfortable wearing sensation. Easy to put on by grip loops. Anatomical knit and closure form ensure precise and comfortable fit, especially in the hip-bending region. Height: ventral 17 cm, dorsal 32 cm. Colors: beige, silvergrey.

INDICATIONS

Application mainly for slipped vertebral disk. Post-operative for clearing of lordosis (e.g. after slipped disk operation), facet syndrome, lumbar spinal duct stenosis, symptomatic stenosis of the intervertebral foramina.

Sizes	Circumference of waist
small	up to 89 cm
medium	89 – 104 cm
large	104 – 119 cm
x-large	119 – 134 cm

1 cm = 0,3937 inch




No. 104 100
BORT
 Abdominal Support

Soft lining at bodyside, comfortable wearing features. Individually adjustable with Velcro closure.
 Height: 21 cm.
 Color: beige.



No. 104 150
BORT PostOban®
 Thoracal-Abdominal Support

Tissue air permeable with large Velcro closure.
 Bipartite hook strap for simplified putting-on and improved fit of the bandage.
 Height: 26 cm.
 Color: white.



Measuring point:
 Circumference waist
 above Crista iliaca.

INDICATIONS

After operations of abdomen, peritoneum injuries. Connective tissue insufficiency. Pendulous abdomen.

INDICATIONS


After abdominal/thoracal operations. Peritoneum injuries. Connective tissue insufficiency. Pendulous abdomen.

Sizes	Circumference of waist
small	up to 87 cm
medium	87 – 102 cm
large	102 – 117 cm
x-large	117 – 132 cm

1 cm = 0,3937 inch

Sizes	Circumference of waist
x-small	65 – 75 cm
small	75 – 85 cm
medium	85 – 95 cm
large	95 – 105 cm
x-large	105 – 115 cm
xx-large	115 – 125 cm

1 cm = 0,3937 inch



No. 104 050
BORT
Scar Rupture Bandage

Elastic, skin-friendly fabric permeable to air with 2 supporting splints. Plastic material reinforced silicone pelotte with adhesive layer for individual positioning.

Height: 13 cm.

Color: white.



Shape of pelotte.



No. 104 070
BORT
Umbilical Hernia Bandage

Distinct shaped silicone pelotte reinforced by plastic material, freely adjustable by adhesive coating.

Elastic, skin-friendly fabric, 2 supporting splints in dorsal area.

Height: 13 cm.

Color: white.



Shape of pelotte.

INDICATIONS

Scar rupture, scar insufficiency, abdominal muscle infirmity/-paralysis.

Sizes	Circumference of waist
-------	------------------------

Size 1	.60 – 95 cm
Size 2	.95 – 125 cm

1 cm = 0,3937 inch


INDICATIONS

Reducible abdominal hernia, as e.g. hernia umbilicalis/paraumbilicalis, hernia ventralis, hernia ventralis lateralis, hernia epigastrica.


Sizes	Circumference of waist
-------	------------------------

Size 1	.60 – 95 cm
Size 2	.95 – 125 cm

1 cm = 0,3937 inch




Measuring point:
Circumference waist
above *Crista iliaca*.




No. 104 620
**BORT Abdominal Support
 for Pregnant Women**

Air-permeable elastic fabric, wide Velcro closure. Additional elastic strap allows adjustment to varying abdomen circumferences. Dorsal relief.
 Height: 24 cm.
 Color: white.




No. 104 090
**BORT
 Stoma Bandage**

Circular abdominal bandage with specific variable stoma recess. Velour laminated material, in dorsal part elastic insert with 2 support splints. Large surface, skin-friendly lining in stoma area to be freely positioned by Velcro fixing.
 Thanks to slashes and Velcro fixing suitable for almost any stoma.
 Washable at 90° C.
 Height: 22 cm.
 Color: white.



Measuring point:
 Circumference pelvis
 above upper Os ilium.



Measuring point:
 Circumference waist
 above Crista iliaca.

INDICATIONS

Rectus diastasis (abdominal wall hernia, connective tissue infirmity), adipositas, multiple pregnancy.

INDICATIONS

Condition after ileostomy and colostomy with concurrent connective tissue insufficiency.

Sizes	Circumference of pelvis
Size 1	90 – 105 cm
Size 2	105 – 120 cm
Size 3	120 – 135 cm

1 cm = 0,3937 inch

Sizes	Circumference of waist
Size 1	70 – 105 cm
Size 2	105 – 135 cm

1 cm = 0,3937 inch




No. 104 160
BORT
 Abdominal Support for Women

Unelastic drill fabric and supporting splints in abdominal and back section.

Height: 34 cm.

Color: beige.



No. 104 170
BORT
 Abdominal Support for Men

Unelastic drill fabric at abdominal and back section. Lateral elastic inserts. 3 supporting splints and buckle closure in abdominal section, cross strap closure.

Height: 27 cm.


Color: beige.

INDICATIONS

After abdomen operations, abdominal wall paresis, abdominal wall hernia.

Sizes	Circumference of hip
Size 1	85 – 90 cm
Size 2	90 – 95 cm
Size 3	95 – 100 cm
Size 4	100 – 105 cm
Size 5	105 – 110 cm
Size 6	110 – 115 cm
Size 7	115 – 120 cm
Size 8	120 – 125 cm
Size 9	125 – 130 cm

1 cm = 0,3937 inch



Measuring point:
 Circumference hip
 above upper front
 Os ilium.

INDICATIONS

After abdomen operations, abdominal wall paresis, abdominal wall hernia.

Sizes	Circumference of hip
Size 1	85 – 90 cm
Size 2	90 – 95 cm
Size 3	95 – 100 cm
Size 4	100 – 105 cm
Size 5	105 – 110 cm
Size 6	110 – 115 cm
Size 7	115 – 120 cm

1 cm = 0,3937 inch



No. 102 700
BORT
Symphysis Belt

Firm, non-elastic pelvis surrounding belt. Inside padded. Conical shape for good fit. Individual and safe fitting by large size Velcro closure. Height: 17 cm. Color: beige.




No. 201 100
BORT StabiloHip
Hip Protector

High quality function panty made of cotton and Lycra. Highly elastic waistband and leg border assure high wearing comfort. Safe positioning of protectors in integrated protector pockets. Easy to take out. Protectors made of visco-elastic soft foam. Special shock absorbing effect. Color: white.

Protector Set (No. 201 120) and Protector Panty (No. 201 200) also available separately.

Measuring point:
 Circumference hip
 above upper front
 Os ilium.



INDICATIONS

Structural disturbance in the area of sacroiliac joint/symphysis pubica, symphysis stretching (symphysiolysis, symphysis dehiscence), symphysis rupture, dysplastic flat joint socket (dysplasia coxae valga, dysplasia coxae congenita).

Sizes	Circumference of hip
small	75 – 90 cm
medium	90 – 110 cm
large	110 – 125 cm

1 cm = 0,3937 inch

INDICATIONS

Tumble of various geneses.

Sizes	Clothing size (D)	
	ladies	men
small	38	46
medium	40	48
large	42	50
x-large	44	52
xx-large	46/48	54/56



No. 105 200
BORT
Hernia Panty for Men

Slip form, with hingeable front part, fully elastic. Can be comfortably used as underpants.
Color: beige.



No. 104 800
BORT
Hernia Aid


Double-sided, pelotte detachable, therefore also single-side usable. With thigh strap.
Color: white.

INDICATIONS

Reducible inguinal hernia, after hernia operation.

Sizes	Circumference of hip
small	80 – 85 cm
medium.	85 – 90 cm
large.	90 – 95 cm
x-large	95 – 100 cm
xx-large	100 – 105 cm

1 cm = 0,3937 inch



Measuring point:
Circumference hip
above upper front
Os ilium.

INDICATIONS

For reducible inguinal hernia.

Sizes	Circumference of hip
small	up to 89 cm
medium.	89 – 104 cm
large.	104 – 119 cm

1 cm = 0,3937 inch



No. 109 100
BORT HerniaFix®
Spring Truss with
Round Pelotte

Spring Truss with inguinal round pressure pelotte. Nonelastic strap with 50 cm long spring. With fabric cover. Single-sided.
Color: beige.



No. 109 200
BORT HerniaFix®
Spring Truss with
Anatomic Pelotte

Spring truss with anatomically shaped pelotte. Nonelastic strap with 50 cm long spring. With thigh strap. Single-sided.
Color: beige.



No. 109 400
BORT HerniaFix®
Elastic Spring Truss with
Round Pelotte, Short Spring

Elastic spring truss with approx. 20 cm **short** spring and round pelotte, single-sided.
Color: beige.

INDICATIONS

Reducible inguinal hernia.



No. 109 700
BORT HerniaFix®
Elastic Spring Truss with Ana-
tomic Pelotte, Short Spring

Elastic spring truss with approx. 20 cm **short** spring and anatomic pelotte, thigh belt, single-sided.
Color: beige.

INDICATIONS

Reducible inguinal hernia.

Sizes	Circumference of pelvis
.....	.85 cm
.....	.90 cm
.....	.95 cm
.....	.100 cm
.....	.105 cm

Please indicate left or right.

1 cm = 0,3937 inch

Sizes	Circumference of pelvis
.....	.85 cm
.....	.90 cm
.....	.95 cm
.....	.100 cm
.....	.105 cm

Please indicate left or right.

1 cm = 0,3937 inch



No. 109 500
BORT HerniaFix®
Elastic Spring Truss,
double-sided

Elastic spring truss with two springs, approx. 20 cm long and 2 round pelottes, double-sided.
Color: beige.



No. 109 800
BORT HerniaFix®
Elastic Spring Truss with Ana-
tomic Pelotte, double-sided


Elastic spring truss with 2 springs, approx. 20 cm long and 2 anatomic pelottes, thigh belt, double-sided.
Color: beige.

INDICATIONS

Reducible inguinal hernia.

Sizes	Circumference of pelvis
.....	.85 cm
.....	.90 cm
.....	.95 cm
.....	.100 cm
.....	.105 cm

1 cm = 0,3937 inch



Measuring point:
Circumference pelvis
above upper Os ilium.

INDICATIONS

Reducible inguinal hernia.

Sizes	Circumference of pelvis
.....	.85 cm
.....	.90 cm
.....	.95 cm
.....	.100 cm
.....	.105 cm

1 cm = 0,3937 inch




No. 110 100
BORT HerniaFix®
Special Truss with Adjustable
Round Pelotte

Special truss with **one** adjustable and exchangeable round pelotte, rubber strap, in the back free of pressure and springs, **single-sided**. Color: beige.



No. 110 300
BORT HerniaFix®
Special Truss with Adjustable
Anatomic Pelotte

Special truss with **one** adjustable and exchangeable anatomic pelotte, rubber strap, in the back free of pressure and springs, **single-sided**, with **thigh belt**. Color: beige.



Measuring point:
Circumference pelvis
above upper Os ilium.

INDICATIONS

Reducible inguinal hernia.

Sizes	Circumference of pelvis
.....	.85 – 94 cm
.....	.95 – 105 cm
.....	106 – 115 cm

Please indicate left or right.

1 cm = 0,3937 inch

INDICATIONS

Reducible inguinal hernia.

Sizes	Circumference of pelvis
.....	.85 – 94 cm
.....	.95 – 105 cm
.....	106 – 115 cm

Please indicate left or right.


1 cm = 0,3937 inch



No. 116 150
BORT
Elastic Thigh Support

Conical shape provides good fit, even for distinct thigh. Flat knitting. Compression class II. 21 cm long. Color: beige.

 **Microfibre comfort**



No. 104 400
BORT
Knee Bandage

Anatomic shape-knitted, therefore solid support and good fit. Corresponds to compression class II. Color: beige.

 **Microfibre comfort**

INDICATIONS

After injuries of the thigh. For treatment after sports injuries.

Sizes	Circumference of thigh
small	up to 43 cm
medium	43 – 50 cm
large	50 – 56 cm
x-large	56 – 61 cm

1 cm = 0,3937 inch


INDICATIONS

For sports injuries, damaged knee ligaments, for articular effusions and **after removal of plaster.**

Sizes	Circumference of knee joint
small	up to 32 cm
medium	32 – 37 cm
large	37 – 42 cm
x-large	42 – 47 cm

1 cm = 0,3937 inch

Measuring point:
e = Circumference middle of knee joint measured over patella.





No. 054 400
BORT
 Dual-Tension Knee Support

Dual-tension elastic, does not slip, thin material, comfortable in the hollow of the knee.
 Compression class II.
 Color: beige.

 **Microfibre comfort**



No. 103 400
BORT
 Knee Support Cross Strap

Crossed straps do not irritate in the hollow of the knee and provide safe support.
 Compression class II.
 Color: beige.

INDICATIONS

For sports injuries, damaged knee ligaments, for articular effusions and **after removal of plaster.**

Sizes	Circumference of knee joint
small	up to 31 cm
medium	31 – 37 cm
large	37 – 42 cm
x-large	42 – 48 cm


1 cm = 0,3937 inch

INDICATIONS


For sports injuries, damaged knee ligaments, for articular effusions and after removal of plaster.

Sizes	Circumference of knee joint
small	up to 32 cm
medium	32 – 37 cm
large	37 – 42 cm
x-large	42 – 47 cm

1 cm = 0,3937 inch



Measuring point:
e = Circumference middle of knee joint measured over patella.



No. 114 410
BORT Stabilo®
 Knee Brace without Stays


Round-knit brace, longwise and crosswise elastic, without lateral stays, compression class II. Patella relief. 30 cm long. Color: beige.

 **Microfibre comfort**

No. 114 450
BORT Stabilo®
 Knee Brace with 2 Stays

Model as above, however with 2 lateral stays each side.

Measuring point:
e = Circumference middle of knee joint measured over patella.



INDICATIONS

For sports injuries, damaged knee ligaments, distortions, after removal of plasters. For arthrosis, meniscus injuries, **damaged lateral ligaments**.

Sizes	Circumference of knee joint
small	up to 34 cm
medium.	34 – 37 cm
large.	37 – 40 cm
x-large	40 – 43 cm


1 cm = 0,3937 inch



No. 210 170
BORT Fillawant
 with Patella Ring

Patella ring of silicone, V-shaped stabilization strips. Thigh parts pulled especially high up, improves muscle function. High portion of cotton.

Recess hollow of knee. Slip-safe at upper and lower edge, versions available with stabilization strips in red (Art. No. 210 170) or in blue (Art. No. 210 160). For definition of size see measuring points (picture).




INDICATIONS

Stabilization of knee joint after surgeries, in case of degenerative signs at ligaments and degenerative changes, rheumatic diseases, arthritis, arthrosis. Chondropathia patellae.

Sizes	Circumference 21 cm above knee	Circumference 15 cm below knee
Size 1	A 39 – 45 B 28 – 33
Size 2	A 45 – 49 B 33 – 36
Size 3	A 49 – 52 B 36 – 39
Size 4	A 52 – 55 B 39 – 42
Size 5	A 55 – 58 B 42 – 45

1 cm = 0,3937 inch




No. 114 150
BORT Knee Bandage with Recess for Patella

- ▶ Patella ring of **silicone**.
- ▶ 1 splint on each side.
- ▶ With pressure-reducing edges.
- ▶ Strong dual-tension quality.
- ▶ Colors: beige, silvergrey.

 **Microfibre comfort**

Measuring points:
 $d = 15\text{ cm}$ below
 $f = 15\text{ cm}$ above
patella.



INDICATIONS

Chondropathia patellae, **instabilities of knee ligaments of different nature**, meniscopathia, post-operative. Bursitis.

Sizes	d-measure	f-measure
small	26 – 29 cm.	32 – 36 cm
medium	29 – 33 cm.	36 – 40 cm
medium plus	29 – 33 cm.	40 – 44 cm
large	33 – 37 cm.	40 – 44 cm
large plus	33 – 37 cm.	44 – 48 cm
x-large	37 – 41 cm.	44 – 48 cm
x-large plus	37 – 41 cm.	48 – 52 cm
xx-large	41 – 45 cm.	52 – 56 cm
xx-large plus	41 – 45 cm.	56 – 60 cm

1 cm = 0,3937 inch



No. 114 420
BORT StabiloGen®


- ▶ Patella fixation bandage.
- ▶ With silicone ring integrated in the lining.
- ▶ Two splints each for improved lateral stability.
- ▶ Strong dual-tension quality.
- ▶ With padding for the hollow of the knee to provide optimal wearing comfort.
- ▶ Colors: beige, silvergrey, blue.

 **Microfibre comfort**

No. 114 520
BORT StabiloGen® Eco

- ▶ Patella fixation bandage.
- ▶ With silicone ring.
- ▶ **One splint each side.**
- ▶ Colors: beige, silvergrey, blue.

Measuring points:
 $d = 15\text{ cm}$ below
 $f = 15\text{ cm}$ above
patella.



INDICATIONS

Patella instability (luxation tendency), articular effusion and swellings, arthrosis and arthritis, post-operative irritations, ligament instability.

Sizes	d-measure	f-measure
For children (No. 114 520, blue only)	23 – 26 cm	28 – 32 cm
small	26 – 29 cm	32 – 36 cm
medium	29 – 33 cm	36 – 40 cm
medium plus	29 – 33 cm	40 – 44 cm
large	33 – 37 cm	40 – 44 cm
large plus	33 – 37 cm	44 – 48 cm
x-large	37 – 41 cm	44 – 48 cm
x-large plus	37 – 41 cm	48 – 52 cm
xx-large	41 – 45 cm	52 – 56 cm
xx-large plus	41 – 45 cm	56 – 60 cm

1 cm = 0,3937 inch




No. 114 390
BORT GenuZip®

- ▶ Patella fixation bandage.
- ▶ With **lateral zip fastener** for easy putting on and taking off.
- ▶ Taking-off shoes is not required to put-on and remove the bandage.
- ▶ Without inconvenient bending down.
- ▶ Integrated silicone ring.
- ▶ One splint each for improved lateral stability.
- ▶ With padding for the hollow of the knee to provide optimal wearing comfort.
- ▶ Strong dual-tension quality.
- ▶ Colors: beige, silvergrey, blue.

 **Microfibre comfort**

*Measuring points:
d = 15 cm below
f = 15 cm above
patella.*




INDICATIONS

Patella instability (luxation tendency), articular effusion and swellings, arthrosis and arthritis, post-operative irritations, ligament instability.

Sizes	d-measure	f-measure
small	26 – 29 cm.	32 – 36 cm
medium	29 – 33 cm.	36 – 40 cm
medium plus	29 – 33 cm.	40 – 44 cm
large	33 – 37 cm.	40 – 44 cm
large plus	33 – 37 cm.	44 – 48 cm
x-large	37 – 41 cm.	44 – 48 cm
x-large plus	37 – 41 cm.	48 – 52 cm
xx-large	41 – 45 cm.	52 – 56 cm
xx-large plus	41 – 45 cm.	56 – 60 cm

1 cm = 0,3937 inch




No. 114 700

BORT Asymmetric[®] based on Dr. R. Weiß

- ▶ Individual guiding of patella towards medial and distal by elastic reins.
- ▶ Anatomically form-knitted.
- ▶ Patella silicone ring.
- ▶ One splint each side.
- ▶ Pressure reduced zones in the peripheral areas and the hollow of the knee.
- ▶ Color: silvergrey/blue.

Measuring points:
 $d = 15\text{ cm}$ below
 $f = 15\text{ cm}$ above
 patella.



INDICATIONS

Patella laterality, luxation tendency, condition after lateral release operation. Patella instability, meniscopathies, chondropathia patellae.

Sizes	d-measure	f-measure
small	26 – 29 cm.	32 – 36 cm
medium	29 – 33 cm.	36 – 40 cm
medium plus	29 – 33 cm.	40 – 44 cm
large	33 – 37 cm.	40 – 44 cm
large plus	33 – 37 cm.	44 – 48 cm
x-large	37 – 41 cm.	44 – 48 cm
x-large plus	37 – 41 cm.	48 – 52 cm
xx-large	41 – 45 cm.	52 – 56 cm
xx-large plus	41 – 45 cm.	56 – 60 cm

Please indicate left or right

1 cm = 0,3937 inch



No. 114 010
BORT Patella
Tendon Strap

Infrapatellar pad inserted in pocket. Special naturally breathing material. Reversed strap for well-dosed pressure.
Washable.
Color: blue/black.



No. 114 610
BORT
Patella Brace

With **fixed** silicone insert.
Medial or lateral usable.
Color: beige.




Measuring point:
Circumference directly
measured below patella.

INDICATIONS

Chondropathia patellae, **femuro-patellar arthrosis.**

Sizes	Circumference of knee joint below patella
Size 1	up to 34 cm
Size 2	34 – 41 cm
Size 3	41 – 48 cm

1 cm = 0,3937 inch



Measuring point:
e = Circumference middle of
knee joint measured over
patella.

INDICATIONS

Chondropathia patellae.
Patella luxation.

Sizes	Circumference of knee joint
small	up to 34 cm
medium.	34 – 39 cm
large.	39 – 44 cm

1 cm = 0,3937 inch




No. 150 140

BORT Stabilo® Knee Bandage with Articulated Joint

- ▶ Strong dual-tension quality with silicone pad.
- ▶ Anatomically form-knitted.
- ▶ Pressure-reduced areas at edges and at hollow of the knee.
- ▶ No seam in the hollow of the knee.
- ▶ Extension limit at 0°.
- ▶ Circular fixation straps for thigh and shank.
- ▶ Color: silvergrey.

Measuring points:
d = 15 cm below
f = 15 cm above
 patella.

**INDICATIONS**

For knee ligament damage, **lateral ligament instability of knee joint**,
 Genu recurvatum.

Sizes	d-measure	f-measure
For children	23 – 26 cm . . .	28 – 32 cm
small	26 – 29 cm . . .	32 – 36 cm
medium	29 – 33 cm . . .	36 – 40 cm
medium plus	29 – 33 cm . . .	40 – 44 cm
large	33 – 37 cm . . .	40 – 44 cm
large plus	33 – 37 cm . . .	44 – 48 cm
x-large	37 – 41 cm . . .	44 – 48 cm
x-large plus	37 – 41 cm . . .	48 – 52 cm
xx-large	41 – 45 cm . . .	52 – 56 cm
xx-large plus	41 – 45 cm . . .	56 – 60 cm

1 cm = 0,3937 inch




No. 114 460

BORT Knee Bandage with Articulated Joint

- ▶ Strong dual-tension quality with silicone pad and **patella recess**.
- ▶ Anatomically form-knitted.
- ▶ Pressure-reduced areas at edges and at hollow of the knee.
- ▶ No seam in the hollow of the knee.
- ▶ Extension limit at 0°.
- ▶ Circular fixation straps for thigh and shank.
- ▶ Color: silvergrey.

Measuring points:
d = 15 cm below
f = 15 cm above
 patella.

**INDICATIONS**

Lateral ligament instability, chondropathia patellae, Genu recurvatum.

Sizes	d-measure	f-measure
small	26 – 29 cm.	32 – 36 cm
medium	29 – 33 cm.	36 – 40 cm
medium plus	29 – 33 cm.	40 – 44 cm
large	33 – 37 cm.	40 – 44 cm
large plus	33 – 37 cm.	44 – 48 cm
x-large	37 – 41 cm.	44 – 48 cm
x-large plus	37 – 41 cm.	48 – 52 cm
xx-large	41 – 45 cm.	52 – 56 cm
xx-large plus	41 – 45 cm.	56 – 60 cm

1 cm = 0,3937 inch



No. 182 150
BORT StabiloPro
Knee Bandage

Dry wearing feeling due to naturally breathing and humidity draining material. Functional guidance/stabilization of patella by special patella ring. Polycentric aluminium gear segment hinge for additional medial and lateral stabilization. With two circular Velcro straps for optimal support.

Adjustable extension at 0°, 10°, 15°, 20°, 30°, 40° and 45°.

Adjustable flexion at 0°, 45°, 60°, 75° and 90°.

Color: black.




No. 182 300
BORT StabiloPro
Knee bandage open form

Dry wearing feeling due to naturally breathing and humidity draining material. Polycentric aluminium gear segment hinge for additional medial and lateral stabilization. Functional guidance/stabilization of patella by special patella ring. With 4 circular Velcro straps for optimal support. Can be opened completely, lockable with Velcro closure.

Adjustable extension at 0°, 10°, 15°, 20°, 30°, 40° and 45°. Adjustable flexion at 0°, 45°, 60°, 75° and 90°.

Color: black.

Measuring point:
 e = Circumference middle of
knee joint measured over
patella.



INDICATIONS

Lateral ligament instability, patella guiding, gonarthrosis.

Sizes Circumference of knee joint

x-small	up to 31 cm
small	31 – 34 cm
medium	34 – 37 cm
large	37 – 40 cm
x-large	40 – 43 cm
xx-large	43 – 46 cm

1 cm = 0,3937 inch

INDICATIONS

Lateral ligament instability, patella guiding. Post-OP, gonarthrosis.

Sizes Circumference of knee joint

x-small	up to 31 cm
small	31 – 34 cm
medium	34 – 37 cm
large	37 – 40 cm
x-large	40 – 43 cm
xx-large	43 – 46 cm

1 cm = 0,3937 inch



No. 114 480
BORT Stabilo®
 Knee Brace Open Style



No. 114 490
BORT Stabilo® Knee Brace
 special width

Bandage can be opened completely.

Ideal for handicapped persons and for injuries not permitting slipping. Steel hinge for good lateral guidance. 4 Velcro closures provide tight support. Color: beige.


With 2 helical springs each side, patella recess and side pads. Especially suitable **for patients with large thighs**. Wedged rubber piece and 2 Velcro closure straps assure non-slipping fit. Color: beige.

INDIKATIONEN

Instability of ligaments.

Sizes	Circumference of knee joint
small	up to 34 cm
medium.	34 – 37 cm
large.	37 – 40 cm
x-large	40 – 43 cm
xx-large	43 – 46 cm

1 cm = 0,3937 inch



Measuring point:
e = Circumference middle of knee joint measured over patella.

INDIKATIONEN

Patella instability (luxation tendency), articular effusion and swellings, arthrosis and arthritis, post-operative irritations, ligament instability.

Sizes	Circumference of knee joint
Size 1	34 – 37 cm
Size 2	37 – 40 cm
Size 3	40 – 43 cm
Size 4	43 – 46 cm
Size 5	46 – 49 cm
Size 6	49 – 52 cm

1 cm = 0,3937 inch



No. 145 000
BORT Immobilization Splint
with Patella Recess

Inner face lined with terry cloth for good wearing comfort. Consisting of 3 parts, medial and lateral splint pocket can be freely positioned. With alu-splints, removable for washing. Velcro closure. Immobilization position in 0° or 20° flexion. Splint length in children's size: 38 cm. Color: black.



No. 145 070
BORT Immobilization Splint
one-piece

One-piece, thus easy to put on. Lateral splint pockets to be freely positioned. Inner side lined with terry cloth. With alu-splints, removable for washing. 0° or 15° flexion. Length approx. 55 cm. Color: black/blue.

INDICATIONS

For complete immobilization of knee joint.

Splint lengths	Fits height
For children	80 – 160 cm
50 cm	up to 160 cm
55 cm	160 – 180 cm
60 cm	over 180 cm

Please indicate flexion 0° or 20° .

1 cm = 0,3937 inch

INDICATIONS

For complete immobilization of knee joint.

One size only.

Please indicate flexion 0° or 15° .



No. 100 600
BORT OTS Knee Brace

- ▶ Optimal stabilization by 4-point-principle.
- ▶ High frame stability in spite of low weight (approx. 600 g) thanks to material-saving design.
- ▶ Paddings available as optional exchange set (No. 100 800).
- ▶ Polycentric gear segment hinge.
- ▶ Limit keys for extension at 0°, 10°, 20°, 30° and 40° and limit keys for flexion at 0° (180°), 45°, 60°, 75° and 90°.
- ▶ Excellent **forming ability** of the brace frame.

*Measuring point:
Width of knee joint at 30°
flexion in position of femur
condyles.*



INDICATIONS

Conservative care of all isolated or combined instabilities of ligament apparatus (e.g. Genu recurvatum). Post-operative after:

- rupture of front cruciate ligament (ACL)
- rupture of rear cruciate ligament (PCL)
- injuries of inner and outer ligament apparatus.

Sizes	Width of knee joint at 30° flexion in position of femur condyles
small	8,5 – 9,5 cm
medium	9,5 – 10,5 cm
large	10,5 – 11,5 cm
x-large	11,5 – 12,5 cm

Please indicate left or right.

1 cm = 0,3937 inch



No. 116 250
BORT
Calf Support

Anatomically knitted.
With seam, skin-friendly.
Compression class II.
Length 23 cm.
Color: beige.


 **Microfibre comfort**



No. 104 500
BORT
Ankle Support

Slip-on, without toe or heel.
Compression class II.
Color: beige.


Measuring point:
Circumference middle
of calf.



INDICATIONS

After injuries of the calf. For
treatment after sports injuries.

Measuring point:
 b = Circumference of ankle
corresp. b measurement
compr. stockings.



INDICATIONS


For injuries of the ankle joint as well
as after operations of the ankle and
of the ankle joint, **for loose liga-
ments in the ankle joint**, subluxa-
tio supinatoria, **pes valgus trouble**
as well as after removal of plaster,
sinew irritations.

Sizes	Circumference of calf
small	up to 33 cm
medium	33 – 38 cm
large	38 – 43 cm
x-large	43 – 47 cm

1 cm = 0,3937 inch

Sizes	Circumference of ankle
small	up to 21 cm
medium	21 – 23 cm
large	23 – 25 cm
x-large	25 – 27 cm

1 cm = 0,3937 inch



No. 054 500
BORT Dual-Tension
Ankle Support

Anatomic shape. Compression class II. Closed heel. Circular knitted, thus without seam.
Color: beige.

 **Microfibre comfort**




No. 182 200
BORT StabiloPren
Ankle Support

Neopren bandage with closed heel. With 2 Velcro closures for variable compression and support of supination or pronation. Limitation of talus advance.

Color: black.

Measuring point:
 b = Circumference of ankle corresp. b measurement comp. stockings.



INDICATIONS

For ankle joint injuries, after ankle and ankle joint operations, subluxatio supinatoria, **pes valgus troubles**, after removal of plaster, **loose ligaments**, sinew irritations.

Sizes	Circumference of ankle
small	up to 21 cm
medium	21 – 23 cm
large	23 – 25 cm
x-large	25 – 27 cm

1 cm = 0,3937 inch

INDICATIONS

Instability of ankle joint area.

Sizes	Shoe sizes (D)
x-small	up to 38
small	38 – 39
medium	39 – 42
large	42 – 46
x-large	46 – 48



No. 054 700
BORT TaloStabil®
Ankle Support
with Silicone Pads

With 2 **structured** silicone pads for intermittent massage. Anatomical shape. Flat knitted. Compression class II. Relief of pressure at ankle ending. Extension at side of foot sole avoids curling.
Colors: beige, silvergrey.



No. 054 600
BORT TaloStabil® Eco
Ankle Support

With 2 silicone pads. Anatomical shape. Flat knitted. Compression class II. Relief of pressure at ankle ending. Extension at side of foot sole avoids curling.
For left and right application.
Colors: beige, blue.

INDICATIONS

Distortions, arthrosis, achillodynia. Ligament instability. After removal of plaster.

Sizes	Circumference of ankle
small	up to 21 cm
medium.	21 – 23 cm
large.	23 – 25 cm
x-large	25 – 27 cm

Please indicate left or right.

1 cm = 0,3937 inch

INDICATIONS


Distortions, arthrosis, achillodynia. Ligament instability, after removal of plaster.

Sizes	Circumference of ankle
small	up to 21 cm
medium.	21 – 23 cm
large.	23 – 25 cm
x-large	25 – 27 cm

Same for left and right.

1 cm = 0,3937 inch





No. 054 900
BORT AchilloStabil®
Ankle Support


Three-dimensional form-knitted dual-tension elastic bandage.
Compression class II.
With especially padded seam.
Viscoelastic profile pads at para-achilles area. Pressure relief of achilles tendon and intermittent massage effect for better blood circulation.
Pressure reduced edge. Additional relief of achilles tendon by special heel wedge (optional).
With heel cushion for non-effected leg to equalize length differences.
Colors: beige, silvergrey.



No. 114 100
BORT Stabilo®
Ankle Support

U-Shaped stabilization element made of plastic material and Stabilo® Velcro strap enable limitation of talus prolapse. Comfortable and dry wearing feeling due to naturally breathing and humidity draining material. Also wearable in sports and walking shoes. Individual adjustment by lacing. Washable.
Colors: white, blue/black.

Measuring point:
 b = Circumference of ankle
corresp. b measurement
compr. stockings.



INDICATIONS

Acute and chronic achillodynia, post-operative (follow-up treatment) after ruptures.

Sizes	Circumference of ankle
small	up to 21 cm
medium	21 – 23 cm
large	23 – 25 cm
x-large	25 – 27 cm

1 cm = 0,3937 inch

INDICATIONS

To stabilize the ankle joint after operations or removal of plaster.

For immobilization in the case of inflammatory irritations, rheumatism, arthritis, overstrained pes valgus. Tendinitides.

Sizes	Circumference of ankle
x-small	up to 19 cm
small	19 – 21 cm
medium	21 – 23 cm
large	23 – 25 cm
x-large	25 – 27 cm

1 cm = 0,3937 inch



No. 100 500
BORT MalleoStabil®-
Orthosis

With two cushion versions:

- Air-/Gel cushion can be pre-cooled. For application during the **acute phase**.
- Moulding foam cushion for optimal wearing comfort during the **rehabilitation phase**.

Anatomically shaped, 25 cm long plastic shells. Distinct stabilization of the upper ankle joint. No impairment of walking.

Color: black.

Also available only with moulding foam cushion (No. 100 520).

INDICATIONS

For stabilization of ankle joint.

For reduced plaster phase. Ligament injuries, upper ankle joint distortions, supination- and/or peripheral rotation traumata. Fibula ligament rupture.

One size only.

Same for left and right.




No. 054 800
BORT TaloFX

Soft-Orthosis made of elastic, humidity draining material. Bilateral shells, soft padded with memory foam to provide high wearing comfort and optimal stability. Positioning of ankle joint in pronation or supination. Limitation of talus prolapsus by cross reins.

With DynamicForceControl (DFC) technology for effective stabilization. Easy to put-on and short design, thus no problem to wear in sports or walking shoes.

Color: black.

Measuring point:
 b = Circumference of ankle
corresp. b measurement
compr. stockings.



INDICATIONS

Conservative treatment of pronation and supination traumas, post-operative protection after fibular ligament rupture, as permanent stabilization aid in case of chronic instability of the upper and/or lower ankle joint.

Sizes

Circumference of ankle

For children	17 – 19 cm
small	19 – 21 cm
medium	21 – 23 cm
large	23 – 25 cm

Please indicate left or right.

1 cm = 0,3937 inch



No. 112 070
BORT Metatarsal
Bandage

Elastic bandage. Pelotte inserted in aperture.

Packing unit: pair.

Also available without pelotte (No. 111 070).

Color: beige.



No. 112 090
BORT Hammer Toe
Correction Bandage

Individually adjustable, light, easy to put on. One bandage is necessary per correction.

Several bandages can be worn at the same time.

Color: white.



Measuring point:

a = Circumference ossa metatarsalia corpus I-V.

INDICATIONS

For flat foot and splayfoot.

Sizes	Circumference cm around forefoot
.....	16, 17, 18, 19, 20,
.....	21, 22, 23, 24, 25, 26

..... 16, 17, 18, 19, 20,
..... 21, 22, 23, 24, 25, 26

1 cm = 0,3937 inch

INDICATIONS

To correct hammer toes.

One size only.



No. 930 010

BORT Valco® Hallux Valgus Bandage

- ▶ The **anatomically adapted shape** of the splint guarantees the physiological correction of the basic joint of the big toe.
- ▶ **Good wearing comfort** thanks to the generously dimensioned foam paddings which safely protect in the critical area of the toe semi-shell.
- ▶ To avoid excessive heat, the splint is provided with **a ventilation hole in the toe area.**
- ▶ **Two adjustable Velcro closures** allow a precise adjustment of the moments of leverage with the splint remaining in the **same optimum position.**
- ▶ Packing unit: **1 piece.**
- ▶ Color: white/blue.

INDIKATIONEN

For hallux valgus, post-operative capsule relief.

Sizes	Shoe sizes (D)
small	up to 38
medium.	38 – 41
large.	over 41

Please indicate left or right.



No. 470 620
BORT
Splint Peroneal

Lightweight orthosis made of polypropylene, also fits into shoe. Highly flexible half-shell in calf area with Velcro closure for individual adjustment to different calf circumferences.

3/4-sole up to the metacarpophalangeal joint of toes.

Color: white.



No. 930 120
BORT
Hallux Valgus Shoe

To be totally opened, extremely soft forefoot, roll-on sole, forefoot relief by negative heel.

Packing unit: **1 piece.**

Color: white.

INDIKATIONEN

Peroneal paralysis.

Sizes	Shoe sizes (D)
small	36 – 39
medium	40 – 43
large	44 – 48


Please indicate left or right.

INDIKATIONEN

After hallux valgus OP, operations at forefoot, as bandage shoe.

	Shoe sizes (D)
.....	35/36
.....	37/38
.....	39/40
.....	41/42
.....	43/44
.....	45/46


Please indicate left or right.



No. 950 100
BORT
Insole

Foot rest area with special foam for optimal adjustment to foot shape. Lower side throughout of visco-elastic material, self-adhesive. Usable as insole with excellent shock absorption or as support for various modules.

Packing unit: **pair**.



No. 950 200
BORT
Heel Cushion

With visco-elastic lower side, self-adhesive, thus ideal for sandals. Foot side textile covered.

Packing unit: **pair**.

And these are the decisive advantages of the BORT insert systems:

- The modules can be glued into the shoes.
- To wear also in open shoes.
- Washable and therefore very hygienic.

INDICATIONS


Shock absorption, damping.

Sizes	Shoe sizes (D)
A	32 – 35
B	36 – 40
C	41 – 44
D	45 – 48

INDICATIONS


Shock absorption, damping.

Sizes	Shoe sizes (D)
A	32 – 35
B	36 – 40
C	41 – 44
D	45 – 48



No. 950 220
BORT
Heel Spur Cushion

With removable cushion module, self-adhesive, thus ideal for sandals. Foot side textile covered.
Packing unit: **pair**.



No. 950 240
BORT
Heel Lift

For compensation of leg length differences. Self-adhesive.
Packing unit: **piece**.

And these are the decisive advantages of the BORT insert systems:

- The modules can be glued into the shoes.
- To wear also in open shoes.
- Washable and therefore very hygienic.

INDICATIONS

Heel spur, shock absorption.

Sizes	Shoe sizes (D)
A	32 – 35
B	36 – 40
C	41 – 44
D	45 – 48

INDICATIONS

For compensation of leg length differences.

Sizes	Shoe sizes (D)
small	up to 40
large	over 40
Heights	5,0 mm or 10,0 mm



No. 930 050
BORT Silicone
Heel Cushion

Heel cushion made of silicone.
Reduces shock load onto ankle,
knee and hip joint as well as onto
spine.
With laterally raised edges.
Height: 6 mm.
Packing unit: **pair**.



No. 930 070
BORT Silicone Heel Spur
Cushion with SoftSpot

Heel cushion made of silicone,
with SoftSpot insert.
Anatomically designed soft bedding
in the heel spur area. Defined and
efficient pressure relief by SoftSpot-
technology. Effective reduction of
shock impact onto ankle, knee and
hip joints, as well as onto spine.
With laterally raised edges.
Height: 6 mm.
Packing unit: **pair**.

INDICATIONS

Shock absorption, damping.

Sizes	Shoe sizes (D)
small	up to 38
medium.	39 – 42
large.	43 – 46

INDICATIONS

Heel spur, shock absorption.

Sizes	Shoe sizes (D)
small	up to 38
medium.	39 – 42
large.	43 – 46


PEDI SOFT®



**TEXTILES
VERTRAUEN**
Schadstoffgeprüfte Textilien
certification no. 03.0.9028 R. Hohenstein
(recognized ecological test procedure)

No. 950 120
PediSoft® Air Cushion Insole

- ▶ Insole comprises a total of 4 in pairs arranged air chambers.
- ▶ By defined air exchange between the individual chambers the naturally rolling action of the foot will be optimally supported.
- ▶ The insole has an insulated to cold effect and a wear-resistant textile surface.
- ▶ Can be washed in the washing machine at a temperature of 30°C.
- ▶ Colors: beige, blue.
- ▶ Packing unit: pair.



Dimensions of display:
height 42,0 cm
width 22,5 cm
depth 19,0 cm



INDICATIONS

Shock absorption, damping. Improvement of sensorimotor function, stimulation of foot reflex zones, reduction of muscular dysbalances.

	Shoe sizes (D)
.....	36/37
.....	38/39
.....	40/41
.....	42/43
.....	44/45



No. 190 500

MultiCure® "S" Anti-Decubitus Gel Cushion

- ▶ Novel solid gel with elasticity, viscosity and flexibility so far not reached for gel types.
- ▶ 600% elasticity for optimal adaption.
- ▶ Not odour absorbing.
- ▶ Machine washable.
- ▶ Sterilizable.
- ▶ Totally non-toxic.
- ▶ High pressure tested.

No. 190 510

**MultiTask®
Anti-Decubitus Gel Cushion
Size 40 x 40 x 2,0 cm**

As art. no. 190 500, however with removable cover.

INDICATIONS

Decubitus prophylaxis and therapy support, cold and warmth therapy.

Sizes for No. 190 500

40 x 40 x 3,0 cm	full cushion
40 x 40 x 3,0 cm	ring-shaped cushion
45 x 40 x 3,0 cm	full cushion
45 x 40 x 3,0 cm	ring-shaped cushion
65 x 40 x 3,5 cm	full cushion

1 cm = 0,3937 inch



No. 190 530
BORT
 Codyn® Seat Cushion


Special ring cushion with **gap in the coccyx area**, inflatable, with valve, suitable for any type of seating accommodation.

Dimensions:

approx. 45 cm x 45 cm.

Material: polyurethane.

Color: blue.



No. 198 500
BORT
 Foot-Lift

Lift-off the foot from resting place by thick foam material cuff. Washable, with Velcro closure.

Color: blue.

INDICATIONS

Coccygodynia, condition after coccyx contusion and fracture, perianal relief of hemorrhoids. Postoperational condition after rectum surgery, peripheral abscess.
Condition after episiotomy.

One size only.

INDICATIONS

Decubitus prophylaxis and therapy support.

One size only.



No. 190 300
BORT Anti-Decubitus
Heel Protector

100% synthetics, washable.
Also applicable as knee or elbow
protection. Slip-type.
Packing unit: 1 piece.
Color: white.



No. 190 800
BORT Pads for
Forearm Walking Sticks

Self-adhesive, textile-lined surface,
shape steady.
Dimensions: 11 x 7,5 cm.
Packing unit: pair.
Colors: black, blue.

No. 190 820
BORT Grip-Padding

Inside neoprene, outside textile
lining. Covered Velcro closure.
Dimensions: 10 x 10 cm.
Packing unit: pair.
Color: black.

INDICATIONS

Decubitus prophylaxis and therapy
support.

Sizes	Circumference of ankle
small	up to 21 cm
medium	21 – 25 cm
large	25 – 29 cm
x-large	29 – 33 cm

1 cm = 0,3937 inch

INDICATIONS

Pressure irritations by extraordinary
load of hand or elbow area.

One size only.



No. 951 500
BORT Putty

Silicone putty for therapeutic finger exercises in 8 different types of hardness, bacteriostatic. 80 g or 440 g box.



No. 112 150
BORT
VITAL-Back Cushion

Special adjustment technique enables:

- infinity variable height adjustment
- 3 depth positioning for support of lordosis

Therefore individual adaptation.

Relief of lumbar spine. Supports ergonomically correct seating position.

Wear-resistant, washable cover.

Available with and without fixing strap.

Dimensions: 42 x 30 x 5 cm.

Colors: blue, black and brown.

INDICATIONS

For regaining mobility of the hand.

INDICATIONS

Back pains of unknown reason, lumbago, protection against wrong positioning during sitting, stresses of back muscles, intervertebral disc problems.

Hardness	Color
super soft	white
extra soft	light blue
soft	turquoise blue
middle soft	light green
middle	dark blue
middle solid	blue green
solid	dark green
very solid	grey



No. 118 150 / 118 170 / 118 190

BORT EasyLife® Reacher


- ▶ With special rubber lips for gripping of large and small objects (bottles, needles, etc.).
- ▶ With the adjustment positions:
 - short leverage – normal force with unrestricted gripping ability of hand.
 - long leverage – therefore lower force - with restricted gripping ability of hand.
- ▶ With locking device, i. e. gripped items will be held also after release of grip handle.

INDICATIONS

e.g. paraplegia, mobility restriction.


	Length
.....	50 cm (No. 118 150)
.....	75 cm (No. 118 170)
.....	90 cm (No. 118 190)

1 cm = 0,3937 inch



No. 118 100
BORT EasyLife®
 Reacher, foldable

With special rubber lips for gripping of large and small objects (bottles, needles, etc.).
 Ergonomically shaped handle.
 With locking device, i. e. gripped items will be held also after release of grip handle.
 Length: 75 cm, folded 40 cm.
 Therefore ideal for anywhere and every situation.



No. 117 650
BORT
 Dressing Help

One end forms a shoe horn, the other end offers a hook with multiple use facilities.
 Length: 60 cm.

INDICATIONS

e.g. paraplegia, mobility restriction.

One length only.

INDICATIONS

Reduced gripping ability, limited mobility.

One length only.



ActiveColor Supports

Joint Supports for Work, Leisure and Sports.



No. 1410
ActiveColor
Wrist Support

Problems at the wrist can occur due to a onetime heavy stress (e.g. straining) as well as due to continued overstrain (inflammation of tendon). The ActiveColor Wrist Support with variable Velcro closure braces and offers individually adjustable compression.

Material: Polyamide, Viscose.


Colors: white (see picture), red, blue, beige, black.

INDICATIONS

Overstrain prophylaxis, sprains and strains, overstress irritations.

Sizes	Circumference of wrist
medium	up to 19 cm
large	over 19 cm

1 cm = 0,3937 inch



No. 1430
ActiveColor
Thigh Support

The ActiveColor Thigh Support prevents overstrain of the thigh muscles. Also for aftertreatment of sports injuries at thigh.

Length: 21 cm.

Material: Polyamide, Lycra.

Colors: beige, blue (see picture).

 **Microfibre comfort**

No. 1420
ActiveColor
Elbow Support

With strong compression force. Recommend for all sports with intensive arm use.

Colors: white, blue, red, black, beige.

INDICATIONS

Overstrain prophylaxis, sprains and strains, overstress irritations.


Sizes No. 1430	Circumference of thigh
small	up to 43 cm
medium	43 – 50 cm
large	50 – 56 cm
x-large	over 56 cm

Sizes No. 1420	Circumference of elbow
small	up to 26 cm
medium	26 – 28 cm
large	over 28 cm


1 cm = 0,3937 inch



With air-permeable storage bag on one's way.



No. 1460
ActiveColor
Knee Support with Zipper



No. 1450
ActiveColor
Ankle Support

With lateral zipper for easier putting on and taking off. One lateral stay each.

Material: Polyamide, Elasthan.
Colors: beige, blue, black (see picture).

Microfibre comfort

Perfect fit by extended pendulum heel. The ActiveColor Ankle Support provides effective bracing and can be worn in sports as well as in walking shoes.

Material: Polyamide, Lycra.
Colors: white, blue, red (see picture), black, beige.

Microfibre comfort

No. 1440
ActiveColor
Knee Support

Without lateral stays and zipper.
Colors: white, red, blue, black, beige.

Dimensions of display:
height 55,0 cm
width 31,0 cm
depth 42,0 cm



INDICATIONS

Overstrain prophylaxis, sprains and strains, overstress irritations.

INDICATIONS

Overstrain prophylaxis, sprains and strains, overstress irritations.

Sizes	Circumference of knee joint
small	up to 32 cm
medium	32 – 37 cm
large	37 – 42 cm
x-large (blue, black, beige)	over 42 cm

1 cm = 0,3937 inch

Sizes	Circumference of ankle
small	up to 21 cm
medium	21 – 23 cm
large	23 – 25 cm
x-large (blue, black, beige)	over 25 cm

1 cm = 0,3937 inch

CLIMA[®] Care

The thermal program with the cosy-warm feeling.

**TEXTILES
VERTRAUEN**
Schadstoffgeprüfte Textilien
Acc. to Öko-Tex standard 100
Certification no. 20.0.4316 R. Hohenstein



No. 106 500
CLIMACare[®]
Shoulder Warmer


Light and inconspicuous, therefore well to wear under the clothing. Protects shoulder and neck areas which are especially endangered during cold and draft.

With it's comfortable thermal effect helps to sooth pains and stress.

Material composition:
80 % cotton, 20 % angora.

Packing unit: **1 piece.**

Color: white.



No. 106 600
CLIMACare[®]
Body Warmer

Very effective protection against cold by high angora portion. For pain in the lumbar vertebra area, for sports and other open-air activities.

Material composition:
44 % angora, 26 % virgin wool,
26 % polyamide, 4 % elasthan.

Packing unit: **1 piece.**

Colors: white, beige.


Sizes	Clothing size (D)	
	ladies	men
small	36/38	48
medium	40/42	50
large	44/46	52
x-large	48/50	54

Sizes	Circumference of hip
small	up to 98 cm
medium	98 – 108 cm
large	108 – 118 cm
x-large (beige only)	118 – 128 cm


1 cm = 0,3937 inch

CLIMA Care

The thermal program
with the cosy-warm
feeling.



No. 106 700
CLIMACare®
Joint Warmer



No. 106 800
CLIMACare®
Foot Warmer

High elastic knitting. Therefore well to wear without restriction of mobility. Size "small" also usable as elbow warmer.

Material composition:

40 % angora, 30 % polyamide,
25 % virgin wool,
5 % Elasthan.

Packing unit: **pair**.

Colors: white, beige.



Dimensions of display:
height 153 cm
width 43 cm
depth 36 cm

Tough and warming virgin wool-
angora mix. Ideal in case of irritation
to fall asleep by cold feet. For more
indoor walking safety by anti-skid
naps.

Material composition:

67 % wool, 25 % angora,
8 % polyamide.

Packing unit: **pair**.

Colors: white, blue.


Sizes	Circumference of joint
small	up to 32 cm
medium	32 – 38 cm
large	38 – 44 cm
x-large (beige only)	44 – 50 cm

1 cm = 0,3937 inch


Sizes	Shoe sizes (D)
x-small (blue only)	31 – 35
small	36 – 38
medium	39 – 40
large	over 41

PEDI SOFT®


Soft Foot Protection made of Silicone. Ideal also for **diabetics**.




No. 107 010
PediSoft®
Toe Spreader
small or large



No. 107 020
PediSoft®
Toe Separator
small or large




No. 107 030
PediSoft®
Bunion Shield




No. 107 060
PediSoft®
Silicone Ring
small or large




No. 107 070
PediSoft®
Heel Cushion
small or large




No. 107 080
PediSoft®
Forefoot Cushion



No. 107 100
PediSoft®
Toe Cushion
right or left



No. 107 110
PediSoft®
Bunion Cap



No. 107 120
PediSoft®
Forefoot Insole

The silicone material offers decisive advantages:

- ▶ It is **soft** and **smooth** just like the own skin and comfortable to wear.
- ▶ Thanks to **memory effect** it adjusts itself optimally to the contour of each part of the foot without slipping and returns to its original shape immediately again after wearing.
- ▶ After use it is washable with soap and can be disinfected if required, and therefore is **more hygienic** as other materials.
- ▶ It is **long-lasting**, wear-resistant and designed for multiple use.



Proven quality
with textile cover.




No. 137 030
PediSoft® TextLine
Hallux-Pad
small or large




No. 137 080
PediSoft® TextLine
Forefoot-Pad
small or large



No. 137 040
PediSoft® TextLine
Forefoot-/Hallux-Pad
small or large




No. 137 060
PediSoft® TextLine
Tube Bandage
small or large



No. 137 050
PediSoft® TextLine
Toe-/Finger Cap
small, medium or large



No. 137 010
PediSoft® TextLine
Toe-/Finger Protection Pad
small, medium or large



No. 137 100
PediSoft® TextLine
Toe-Pad
right or left

Advantages of PediSoft® TextLine:

- ▶ Long lasting, wear-resistant, for multiple use. **Machine washable.**
- ▶ Very hygienic, skin-friendly and comfortable wearing feeling due to **textile cover.**
- ▶ The silicone cushion distributes optimally the pressure and protects against additional external pressure.



Dimensions of display:
height 50,5 cm
width 39,0 cm
depth 20,5 cm

Soft Socks ergo



No. 123 200/123 300/123 400


SoftSocks ergo – Special socks, super soft and compressions-free

- ▶ Excellent quality made of 70% cotton and 30% polyamide.
- ▶ Soft band – no constricting stocking rubber.
- ▶ No noticeable inside seams which could rub.
- ▶ Non-slipping position by special knitting.
- ▶ Extra protection: reinforced heel and closed, manually chain-stitched comfort foot point.
- ▶ High portion of naturally fibre provides for healthy foot climate.
- ▶ Perfect fit by **ergonomically shaped** right and left socks.
- ▶ **Extra-wide version** for ankle circumference up to 33 cm.
- ▶ Also suitable for diabetics and people with vascular disease.
- ▶ Made in Germany.
- ▶ Colors: blue, silvergrey, black, beige.

Dimensions of display:
height 50,0 cm
width 48,0 cm
depth 26,0 cm



Measuring point:
 b = Circumference of
ankle corresp.
 b measurement
compr. stockings



Sizes	Shoe Sizes (D)
.....	35 – 37
.....	38 – 40
.....	41 – 43
.....	44 – 46
Width	Circumference of ankle
normal	up to 23 cm
wide	23 – 27 cm
extra-wide (black and beige only)	28 – 33 cm

1 cm = 0,3937 inch



No. 520 000
StabiloColor®

Available in colors yellow, green, red, blue, violet, as well as in the widths 4 cm, 6 cm, 8 cm, 10 cm. Length: 5 m (stretched).
Avoid the impression of illness and disability by bright colors. Have a positive psychological effect on grown-up and little patients.
They are wash-proof, without losing their shiny colors and elasticity.
Top quality **with woven edges** made of 88% cotton, 11% polyamide and 1% polyurethane.
Elasticity approx. 80%.



No. 521 000
StabiloColor® haft

Cohesive version.
Colors and widths as adjoining.



No. 522 000
StabiloColor® haft
latex-free

Latex-free, cohesive version.
60% cotton, 39% polyamide, 1% polyurethane.
Colors. beige, blue.
Widths: 6 cm, 8 cm.
Length: 5 m stretched.

Dimensions
of display:
height 32,0 cm
width 33,0 cm
depth 39,0 cm



PEDI SOFT® BLISTER PAD

For prophylaxis and therapy of blisters.



No. 107 620
PediSoft®
Blister Pad

Protects like a second skin. Supports the natural healing process of blisters. If necessary the hydrocolloidal wound padding absorbs wound secretion and forms a padding gel core which protects against pressure and friction.

Contents: 5 pieces.
Dimensions: 6 x 4 cm.

One size only.

INDICATIONS

For prevention and treatment of blisters.

Dimensions of display:
height 34,0 cm
width 25,0 cm
depth 24,0 cm



PEDI SOFT® FOOT BALM

Care and protection for sensible skin. Suitable for diabetics



No. 107 410
PediSoft®
Foot Balm

Intensive care for very dry, chapped and sensitive skin. Suitable for diabetics. The active camomile substance has an anti-inflammatory and soothing effect, the meristem extract increases the resistance of the skin, allantoin cares for the skin, making it soft and elastic. Panthenol regenerates, urea absorbs humidity and makes the skin soft. Contents: 50 ml.

INDICATIONS


Prevention and care for dry, chapped and sensitive skin. Especially suitable also for diabetics.

Dimensions of display:
height 57,5 cm
width 15,5 cm
depth 26,0 cm




 EasyLife®

The useful assistants of BORT


 M D M D F S S

No. 108 090
EasyLife®
7-Day-Pill-Box

Clearly arranged 7-Day-Box for controlled and safe taking of medicine and vitamins.

Easy to open, with braille.

Colors: transparent, blue.

Dimensions: 16 x 3 x 2 cm.



No. 108 080
EasyLife®
Medicine Organizer

Clearly arranged 7-Day-Organizer with 4 sections per day. Day compartments separately removable to take along. With braille.

Colors: white, blue.


Dimensions: 17 x 11 x 2 cm.



No. 108 050
EasyLife®
Pill Splitter/Crusher

Easy and safe splitting and pulverizing of pills, medicaments and vitamins. Sharp, stainless-steel blade. Easy to clean.

Color: blue.



No. 108 040
EasyLife®
Pill Splitter

Easy splitting of pills, medicaments, vitamins. Sharp stainless-steel blade. Easy to clean.

Colors: blue, transparent.

Collection

newlook

The reading aid collection
with fashionable flair.

No. 129 100

NewLook Spectra

Diopter +1,0 +1,5 +2,0 +2,5 +3,0



No. 129 900

NewLook Classic


Diopter +1,0 +1,5 +2,0 +2,5 +3,0



No. 129 200

NewLook Crystal

Diopter +1,0 +1,5 +2,0 +2,5 +3,0



No. 130 000

NewLook Matrix

Diopter +1,0 +1,5 +2,0 +2,5 +3,0



No. 129 400

NewLook Avalon

Diopter +1,0 +1,5 +2,0 +2,5 +3,0



No. 129 700

NewLook Sun


Diopter +1,0 +1,5 +2,0 +2,5 +3,0



No. 129 500

NewLook Furore

Diopter +1,0 +1,5 +2,0 +2,5 +3,0



No. 129 800

NewLook Sky


Diopter +1,0 +1,5 +2,0 +2,5 +3,0

Dimensions
of display:
height 44,0 cm
width 29,0 cm
depth 20,0 cm



BORT Foot-Roller

The natural relaxing program.



No. 112 200
BORT VITAL
Foot-Roller

With 51 smooth rollers arranged in a staggered mode.

24 x 22 cm, 6 cm high, 950 g weight.



No. 112 180
BORT ERGO
Foot-Roller


Ergonomically shaped with separate rollers for reflection zone and heel massage, 28 x 21 cm, 8 cm high, 980 g weight.



No. 112 190
BORT
MINI-Roller

With massage grooves and rubber rollers at the ends, glazed.

23 cm long, 5 cm diameter, 165 g.



No. 112 210
BORT STABILO®
Foot-Roller

For foot and hand massage.

With 6 specially structured wood rollers, 32 x 11 cm, 5 cm high, 590 g weight.

BSOS

Bort-Swiss Orthopedic Supply, LLC

2221 Applebrook Dr., Commerce, Mi. 48382

248-242-6907 office/ 248-242-6908 fax

beth.swiss@bort-swissortho.com

www.bort-swissortho.com

www.bort.com



Bort GmbH * P.O. Box 1330 * D-71367 Weinstadt, Germany

Ziegeleistraße 39 - 43 * D-71384 Weinstadt-Benzach

Phone **49-(71 51) 9 92 00-0 * Fax **49-(71 51) 9 92 00-50

medical@bort.com

www.bort.com